



## Poster Number: G27

- Cystic Fibrosis (CF) is caused by mutations in both copies of the *CFTR* gene, and affects ~1 in 3000-4000 US births; more than 2000 *CFTR* variants have been identified, with variable frequency within different population groups.
- We developed a streamlined, two-tube PCR/CE prototype assay which detects 67 variants commonly found in the diverse US population.
- Evaluation of AmpliDeX® PCR/CE *CFTR*\* prototype assay produced consistent results for 73 residual clinical samples and eight controls across two laboratories, using a faster and simpler workflow than many other assays.

Seventy-three previously genotyped residual clinical samples and eight controls (multi-variant mixes of cell lines and synthetic gBlock templates) were assayed by two laboratories. Samples were amplified in a two-tube assay using AmpliDeX® PCR/CE *CFTR* Prototype reagents. Allele-specific PCR generates mobility- and dye-tagged amplicons that are resolved by variant type. PCR products were resolved by capillary electrophoresis (CE) on the Applied Biosystems™ 3500xL; one site also analyzed a subset of samples on 3130xl, 3730xl, and SeqStudio™ Genetic Analyzers. Genotypes were determined by prototype modules in GeneMapper 5.0, and compared using automated peak calling and classification software at one site. Each sample batch was processed within five hours from PCR setup through CE and peak calling.



Sample Agreement with Orthogonal Assay			xTAG® Cystic Fibrosis (CFTR) 60 Kit v2			Overall Sample Agreement
			Homozygous or Compound Heterozygous		Heterozygous	
			WT/WT	MUT/MUT	MUT/WT	
Prototype CFTR PCR/CE Assay	Homozygous or Compound Het	WT/WT	21	0	0	21/21(100%)
		MUT/MUT	0	18	0	18/18 (100%)
	Heterozygous	MUT/WT	0	0	34	34/34 (100%)



- We developed and evaluated an AmpliDeX PCR/CE *CFTR* prototype assay that covers 93% variants from the diverse US population, including all targeted mutations recommended by ACMG/ACOG guidelines.
- The prototype AmpliDeX PCR/CE *CFTR* assay has a rapid workflow, requiring only five hours from sample-to-answer with the flexibility to support high-throughput screening or lower-volume diagnostic applications.
- We observed 100% genotype agreement to reference results from 73 residual blood samples at two labs, including SNPs, INDELs, CNVs, and the IVS8 poly-T/TG modifier.
- Prototype automated software results were 100% concordant with reference results and incorporated multiple QC checks, matched-tube assignment, and automated peak calling/classification.