



Two Complementary and Scalable PCR-based Workflows Enable Next Generation Sequencing of Cancer-Associated Genes in FFPE Tumor DNA

Gary J. Latham, Elizabeth Mambo, Liangjing Chen, Andrew Hadd, Jeffrey Houghton, Kalyan Buddavarapu, Tiffany Sanford, Julie Krosting, Alex Adai, Lana Garmire, and Ashish Choudhary. Asuragen, Inc., Austin, Texas 78744

SUMMARY

- Next generation sequencing (NGS) has the potential to report low abundance, clinically actionable mutations in heterogeneous, real-world tumor specimens.
- Multiplex PCR enrichment methods were developed with both focused (16 genes, 35 amplicon) and broad (52 genes, 981 amplicon) cancer gene panels for direct sequencing on two orthogonal NGS platforms, the Illumina GAIIX and the Ion Torrent Personal Genome Machine (PGM).
- The two enrichment methods enabled uniform read coverage, high depth sequencing (>1000X), high sensitivity mutation detection (1-3%), and excellent concordance with Sanger sequencing and other mutation confirmation methods using FFPE samples from different tumor types and block ages.
- The results support the utility of sensitive, accurate, and high resolution mutation profiling across dozens to thousands of loci in FFPE tumor specimens.

MATERIALS AND METHODS

Two FFPE-compatible PCR-based enrichment panels were developed. The first was a multiplexed PCR assay that targeted 35 amplicons in 16 cancer genes, including the most common mutations in the MAPK/ERK and PI3K/AKT pathways. The second included nearly 1000 amplicons from 52 cancer genes. Primers were designed to avoid known SNPs, repetitive sequences, and pseudogenes whenever possible, and included adaptor sequences to enable direct sequencing on either the Ion Torrent PGM or the Illumina GAIIX. Samples included cancer cell lines (ATCC) and residual FFPE blocks from various sources collected up to 14 years prior and represented colon, melanoma, breast, and other cancers. FFPE DNA inputs ranging from 10 ng to 2 ug were evaluated for PCR enrichment, and samples were barcoded up to 36plex/lane (GAIIX) or 13plex/chip (PGM). Workflows for GAIIX NGS required ~3-6 weeks, whereas the time from purified DNA to processed NGS data on the PGM was <3 days.




Figure 1. Two distinct PCR enrichment workflows with NGS applications ranging from broad content screening to molecular diagnostics.

16 Gene Cancer Panel

Oncogene Panel			
ABL1	FGFR1	HRAS	MET
AKT1	FGFR3	JAK2	NRAS
BRAF	FLT3	KIT	PDGFRA
EGFR	RET	KRAS	PIK3CA

Table 1. A 16 gene, 35 amplicon panel represents >95% of these gene mutations in COSMIC.




Figure 3. Titration of FFPE tumor DNA reveals dose-dependent detection of mutations to 1-2% of total NGS reads.

Illumina GAIIX	Sanger Sequencing	
	POS	Wt
POS	19	3
Wt	2	92

Table 3. Illumina GAIIX NGS of 39 FFPE tumor specimens demonstrates 96% concordance with Sanger sequencing.

DNA Sample	Gene with known mutations	Known Codon Change	Type	Mixing ratio	Expected % Mutation	Ion Torrent PGM % Recovered Coverage	Ion Torrent PGM % Read Coverage	Illumina GAIIX % Recovered Coverage
A549	KRAS	G12S	HOM	35%	35	31.9	971	29.2
MIA PaCa-2	KRAS	G12C	HOM	20%	20	22.6	971	17.8
T24	HRAS	G12V	HOM	10%	10	10.9	248	13.2
RKO	BRAF	V600E	HET	15%	7.5	5.9	801	6.5
SK-Mel-2	PIK3CA	H1047R	HET	15%	10.5	6.6	2290	7.9
SK-Mel-2	NRAS	Q61R	HOM	7%	7	10.7	600	11.2
GP2d	PIK3CA	H1047L	HET	5%	2.5	3.5	2290	6.8
GP2d	KRAS	G12D	HET	5%	2.5	3.7	991	4.0
HCT 116	KRAS	G13D	HET	6%	3	1.5	1024	<0.5
HCT 116	PIK3CA	H1047R	HET	6%	10.5	6.6	2290	7.9
SW1116	KRAS	G12A	HET	2%	1	1.4	991	1.3

Table 4. Ion Torrent PGM and Illumina GAIIX NGS quantitatively recovers known mutations from pooled cancer cell DNA.

Ion Torrent PGM	Orthogonal Confirmation	
	POS	Wt
POS	15	0
Wt	1	91

Table 5. Ion Torrent PGM NGS of 16 FFPE tumor specimens demonstrates 99% concordance with confirmation assays, including Luminex Signature®, Illumina and Sanger sequencing.

Cell Line	Gene with Known Mutation	Known Codon Change	Expected Percent Mutation	NGS Percent Mutation
GDP2	PIK3CA	H1047L	5%	6.1%
HCT116	PIK3CA	H1047R	1.5%	1.4%




Figure 2. As few as 1.5% variants can be accurately quantified by ultra deep NGS using the Illumina GAIIX: co-detection of PIK3CA variants at the same locus.




Table 2. Both intact cell line DNA and FFPE DNA demonstrate A) uniform read coverage, and represent B) low levels of "background" base substitution variants (0.3%).




Figure 4. High depth NGS of BRAF amplicons clearly delineates mutation positive samples, including those that are negative by Sanger sequencing.




Figure 5. Representative variant profile of mutation-positive and negative FFPE tumor DNA revealed by Ion Torrent NGS. Mutations identified on the Illumina GAIIX were confirmed on the Ion Torrent PGM.

52 Gene Cancer Panel

Oncogene and Tumor Suppressor Panel					
ABL1	DNMT3A	GNAQ	MET	PTCH1	TP53
AKT1	EGFR	HIF1A	MPL	PTEN	VHL
AKT2	ERBB2	HRAS	NF2	PTPN11	
BRAF	FES	IDH1	NOTCH1	RB1	
CDH1	FGFR1	IDH2	NPM1	RET	
CDK4	FGFR3	IKBKB	NRAS	SMAD4	
CDKN2A	FLT3	JAK2	PAX5	SMARCB1	
CEBPA	FOXL2	KIT	PDGFRA	SMO	
CREBBP	GATA1	KRAS	PIK3CA	SRC	
CTNNB1	GNA11	MEN1	PIK3R1	STK11	

Figure 6. A 52 gene, 981 amplicon cancer panel compatible with massively parallel picodroplet PCR. Amplicons designed for genes shown in bold were sequenced across all coding exons.

Metric	CRC FFPE1	CRC FFPE1 (replicate)
Total Reads	7,654,272	5,601,478
Median reads/amplicon	6375	4471
2-fold of median	57%	57%
5-fold of median	91%	91%



Metric	CRC FFPE2 (1 ug)	CRC FFPE2 (2 ug)
Total Reads	3,180,016	7,835,776
Median reads/amplicon	2682	6438
2-fold of median	65%	62%
5-fold of median	96%	95%





Figure 7. Picodroplet PCR of FFPE DNA using a 981 amplicon cancer panel supports reproducible, ultra deep NGS using the Illumina GAIIX.

Figure 8. Variant fractions detected by NGS in matched frozen and FFPE tumor DNA are highly correlated. Variants were identified by Illumina GAIIX NGS following RainDance enrichment of 981 amplicons spanning 52 cancer genes.



CONCLUSIONS

- Two distinct PCR workflows enabled high depth and enabled uniform enrichment and high depth NGS of cancer-associated gene regions in FFPE DNA from residual clinical specimens.
- Mutation loads as low as 1-3% were accurately identified in both cancer cell line and FFPE tumor DNA; "background" variant detection was only 0.3%. These results have important implications for detection of low-level mutations, such as drug-resistant mutations, that may be clinically relevant.
- PCR-based enrichment of cancer gene "hotspots" in FFPE tumors revealed 96-99% concordance with orthogonal reference methods, including Sanger sequencing.
- Ion Torrent NGS successfully confirmed novel mutations from screening studies using the Illumina GAIIX, suggesting utility for high sensitivity orthogonal mutation confirmation using distinct enrichment procedures and NGS sequencing chemistries.
- The two proposed NGS approaches can accommodate both large-scale, whole exon mutation assessments for ~96 samples per run, as well as "hotspot" mutation analyses across more than a dozen genes with a rapid turnaround time (<1 week).