A NOVEL RT-ddPCR ASSAY ENABLES SENSITIVE AND ACCURATE QUANTIFICATION OF MET EXON 14 SKIPPING IN LUNG ADENOCARCINOMAS

Michael Dodge¹, Stela Filipovic-Sadic¹, Liangjing Chen¹, Shobha Gokul¹, Junya Fujimoto², Vassiliki A Papadimitrakopoulou³, Ignacio I Wistuba², and Gary J Latham¹ ¹Asuragen, Inc., Austin, TX, ²Department of Translational Molecular Pathology, ³Department of Thoracic/Head and Neck Medical Oncology, The University of Texas MD Anderson Cancer Center, Houston, TX

SUMMARY

- MET experiences exon 14 loss (METex14) in 3-4% of lung adenocarcinomas, thereby promoting protein stabilization and oncogenic pathway activation.
- Exon 14 loss or elevated MET levels are clinically actionable events susceptible to kinase inhibitors, including crizotinib and cabozantinib.
- Droplet Digital[™] PCR (ddPCR[™]) can provide accurate copy number information of both wildtype and METex14 status using a rapid and accessible workflow.
- We describe a prototype multiplexed ddPCR assay that quantifies an internal reference gene and the fraction of METex14 skipping.

INTRODUCTION

The proto-oncogene MET experiences donor/acceptor splice site mutations that can give rise to transcripts lacking exon 14 (METex14)¹. Although this "skipping" event occurs at frequencies comparable to ALK rearrangements in lung adenocarcinomas, current METex14 detection technologies are underdeveloped and burdened by a combination of high costs, low resolution, and/or complex workflows. Herein we report the development of a novel, accurate, and accessible 2-channel Droplet Digital PCR (ddPCR) assay that can quantify the fraction of skipped METex14 transcripts in patient samples. Successful identification of this event correlates with MET inhibitor responsiveness, potentially benefiting treatment decisions and outcomes²

MATERIALS AND METHODS

A METex14 RT-ddPCR assay was developed using optimized reverse transcription and PCR reagents, with an optional pre-amplification step incorporated for samples with low nucleic acid quantity or quality. Gene-specific primers and junction-spanning hydrolysis probes were designed for MET exons 13/14 or 14/15 (wildtype (WT), HEX labeled) and 13/15 (METex14, FAM labeled), respectively. Droplets were generated and run on the QX200™ ddPCR platform (Bio-Rad) with a 4-bin output. Samples consisted of cell lines and a set of 125 lung adenocarcinoma FFPEs (MD Anderson) that were sequenced in parallel for METex14 using the QuantideX® NGS RNA Lung Cancer Kit* (Asuragen)



Figure 1. Assay design for detection of METex14 status. A) Pertinent MET genomic intron/exon junctions. Wide-ranging indels and point mutations cluster around donor & acceptor sites, leading to excision of exon 14 during splicing and the subsequent removal of a negative CBL regulation site codon (Y1003). B) Assay primer & fluorescent probe placement options for detecting processed transcripts.

RESULTS

MET exon 14 skipping status was determined with 100% accuracy in representative cell lines from as little as 2 ng RNA inputs. Known copy number titrations demonstrated analytical sensitivity down to 1% METex14 skipped/WT levels, exceeding the sensitivity typically achieved by other approaches. METex14 positive and negative samples were correctly identified from previously characterized FFPE specimens, including two samples that showed 27% and 95% METex14 skipping ratios, respectively.

2A		Copie	s/20uL		
	Probe	13/15	13/14		Accepted
	Sample	Skipped	WT	% Skipped	Droplets
	H596- 50ng	80,600	350	99.6%	9,996
Positive cell line	H596- 10ng	12,860	54	99.6%	13,203
	H596- 2ng	2,900	14	99.5%	13,700
	RT112-50ng	0	5,180	0.0%	15,465
Wildtype	RT112-10ng	0	1,068	0.0%	12,777
	RT112-2ng	0	184	0.0%	14,243
Cantrala	RT-NTC	0	0	-	13,336
Controls	PCR-NTC	0	0	-	11,131

H596-2ng	3,680	0	100.0%	100%	
H596-2ng	3,700	1	99.9%		
RT112-30ng	1	2,900	0.1%	0.10/	
RT112-30ng	0	3,060	0.0%	0.1%	
25% Positive	1,040	2,480	29.5%	20 / 0/	
25% Positive	858	2,240	27.7%	28.0%	
12% Positive	446	2,500	15.1%	15.4%	
12% Positive	484	2,600	15.7%		
6% Positive	202	2,680	7.0%	6.9%	
6% Positive	226	3,060	6.9%		
1.2% Positive	42	2,900	1.4%	1 4 9/	
1.2% Positive	58	3,140	1.8%	1.0 %	
RT-NTC	0	0	-		
RT-NTC	0	0	-	-	



Figure 2. METex14 ddPCR detects low-frequency skipping events with high specificity. A) Absolute ng input detection for H596 (METex14 positive) and RT112 (wildtype MET) cells. B) METex14 (H596) cDNA was titrated into a wildtype (RT112) background based upon levels established in the originating samples. C) Visualization of select data from B.



Figure 3. An endogenous reference gene can be incorporated into multiplex METex14 ddPCR or qPCR assay configurations. A) Inclusion of a stable endogenous reference gene as an internal control. Two separate wells detect the presence of wildtype(14/15)/skipped(13/15) MET and wildtype/control transcripts, respectively. Shared MET 14/15 quantitation serves as an additional cross-reaction control. **B**) As an alternative to ddPCR, a single-reaction 4-channel qPCR can be performed with ddPCR probe/primer components, including the endogenous control. No off targeting was noted.



A	3' Fusion Genes	#
	ALK	
	ROS1	
	RET	
	FGFR3	
	NTRK3	
	NTRK1	
		_



Figure 5. Accurate detection of METex14 positive and negative samples from residual FFPE lung cancer specimens. A) 125 lung adenocarcinoma FFPE samples (MD Anderson) were screened using the QuantideX RNA Lung Cancer Kit, resulting in the identification of two METex14 positives. Select NGS panel content, including relevant MET transcript isoforms, is shown. METe5:e6 is additionally assessed. B) ddPCR data for representative FFPE negatives (samples ADC3 & ADC4) and the two identified positives (ADC7 & ADC15) following pre-amplification and a 1:100 fold input dilution. Percent skipped vs. wildtype MET noted.

CONCLUSIONS

- Kit (RUO)*.

References

Res. (2006). 66: 283-289.

Figure 4. METex14 ddPCR optimizations increase cDNA inputs into the PCR stage by at least 100-fold. A) Increasing RT loading by as much as 30% of the ddPCR volume was well tolerated. Wildtype MET FFPEs were tested from Asuragen's Tumor Bank (14/15 probe). B) An optional preamplification stage significantly improved sensitivity while preserving allele ratios.

• The prevalence and druggability of METex14 skipping in NSCLC makes this lesion a high-priority diagnostic target for precision medicine.

• We describe a prototype ddPCR assay that can accurately quantify the fraction of METex14 within a single multiplexed reaction.

• Further inclusion of a stable reference gene enables relative quantification of wildtype MET transcript expression independent of mechanism, and can be accommodated by both ddPCR and gPCR formats.

• The assay accepts challenging clinical specimen inputs (e.g. FFPEs) and confirms METex14 status as identified using the Asuragen QuantideX NGS RNA Lung Cancer

1. Kong-Beltran et al. Somatic mutations lead to an oncogenic deletion of MET in lung cancer. Cancer

2. Paik et al. Response to MET inhibitors in patients with stage IV lung adenocarcinomas harboring MET mutations causing exon 14 skipping. Cancer Discov. (2015). 8: 842-849.



