Multi-site Evaluation of the AmplideX[®] PCR/CE TOMM40 Kit (RUO) for Rapid and Accurate Genotyping of Poly-T Length Polymorphisms at rs10524523 of the TOMM40 Gene

Sarah Statt^{1*}, EunRan Suh^{2*}, Julie R Thibert¹, Vivianna M Van Deerlin², and Gary J Latham¹ ¹Asuragen, Inc., Austin, TX; ²Dept. of Pathology and Lab Medicine, Perelman School of Medicine at the University of Pennsylvania, Philadelphia, PA

Summary

- Alzheimer's disease (AD) is the most common cause of dementia, accounting for 60-80% of all dementias and impacting 5.4M Americans.
- An improved risk test may provide insights for AD clinical research, having the potential to improve patient management in the future.
- A poly-T repeat in the TOMM40 gene is associated with age-of-onset risk of AD in some studies. However, this region is very AT-rich and resists efficient, highfidelity PCR amplification and detection.
- The AmplideX® PCR/CE TOMM40 Kit (RUO) achieved rapid and reliable genotyping of the rs10524523 poly-T polymorphism with single-nucleotide resolution across multiple operators, days, instruments and laboratories, and in different sample types.

Introduction

There is increasing interest in identifying Alzheimer's disease (AD) populations that can benefit from early intervention. Poly-T polymorphisms at rs10524523 ('523) of the TOMM40 gene have been reported to influence age of onset in late-onset AD (LOAD)¹ and the rate of cognitive decline². Studies that have examined this locus, however, have been limited by assays with cumbersome workflows and/or suboptimal resolution. Here we describe a two-site evaluation of a fast, simple and accurate PCR assay that reproducibly reports TOMM40 poly-T length.

Materials and Methods

The assay was evaluated at both Asuragen (Site 1) and the University of Pennsylvania (Site 2). Both sites genotyped a common set of cell-line samples. Site 1 assessed additional cell-line materials and 14 matched buccal and peripheral blood samples. Site 2 assessed 93 brain tissue samples, covering 2 cohorts. The AD cohort had 64 neuropathologically confirmed AD cases, of which 11 had matching peripheral blood samples, while the control cohort had 29 without AD neuropathology. Sample gDNA was PCR amplified using the AmplideX® PCR/CE TOMM40 Kit (RUO) and amplicons were resolved by capillary electrophoresis (CE) on the 3500xL instrument (Thermo Fisher Scientific). Genotypes were determined from the mobility of target peaks relative to a calibration curve.



Figure 1. TOMM40 gene and risk alleles. TOMM40 (Translocase Of Outer Mitochondrial Membrane 40) is a protein-coding gene that is a part of the translocase of the outer mitochondrial membrane pore subunit. A variable poly-T tract, rs10524523 ('523), located in intron 6 of the TOMM40 gene (blue star), was found to have a strong relationship with the age of onset, further refining the APOE association with AD1. This poly-1 polymorphism was binned into the following TOMM40 genotypes: Short (S, ≤19 Ts), Long (L, 20-29Ts), and Very Long (VL, ≥30 Ts), with L showing the highest risk¹.



Figure 2 Time-motion for the AmplideX PCR/CE TOMM40 (RIIO) assay workflow An overview of the TOMM40 assay workflow, highlighting the rapid turnaround time of -2.5 hours from gDNA to answer. Total hands on time is ~30 min for 24 samples

Results



Figure 3. The TOMM40 PCR/CE assay produces reproducible poly-T genotypes across multiple operators and days. DNA samples with 16-36 Ts were assessed by three operators across two sites, three different days and three 3500xL instruments. All results were in agreement within a single nucleotide repeat. Representative electropherograms are shown for A) an S/L sample of 16/29 Ts and B) an L/VL sample of 29/36 Ts.



Figure 4. TOMM40 genotypes generated at two different laboratories are identical across 16 cell-line/ whole blood gDNA samples. Evaluation of the AmplideX PCR/CE TOMM40 Kit at both Asuragen (Site 1) and the University of Pennsylvania (Site 2) resulted in 100% agreement for 16 samples. An example of the output from both sites is shown for two blood samples.

	APOE genotype		Age at Onset	Age at Death	TOMM40 Brain	TOMM40 Blood	TOMM40 genotype
Sample1	E3/E3	м	76	85	16/16	16/16	S/S
Sample2	E3/E3	F	83	90	36/37	36/37	VL/VL
Sample3	E3/E3	F	67	78	36/37	36/37	VL/VL
Sample4	E3/E4	м	85	87	16/34	16/34	S/VL
Sample5	E3/E4	м	75	79	16/29	16/29	S/L
Sample6	E3/E4	м	75	86	16/29	16/29	S/L
Sample7	E3/E4	F	66	85	29/34	29/34	L/VL
Sample8	E4/E4	F	N/A	80	29/29	29/29	L/L
Sample9	E4/E4	м	66	72	22/29	22/29	L/L
Sample10	E4/E4	м	74	77	29/29	29/29	L/L
Sample11	E4/E4	F	79	85	29/30	29/30	L/VL
100000444 Assessy 500 110 1000 1000 1000 1000 1000 0 0 0 1000	29T	34T	Brain	TOMARKA JAnnar 50 10 2000 - 2000 - 1000 -	29T	30T Sa	Brain mple11
10000148 Assay 1000 110 2000 1000 1000 1000 1000 0 0	2 [®] T→	3 ⁴ T	Blood	TOOMS-00_Assey	29T	30T Sa	Blood

Figure 5. TOMM40 genotypes agree in matched brain and whole blood samples. Site 2 demonstrated full agreement in the number of poly-T repeats for 11 paired samples, with a median time difference between paired collections of 1520 days. Sample7 and Sample11 electropherograms are shown as an example of the genotype profiles. A separate sample set also revealed genotype concordance for matched buccal and whole blood samples at Site 1 (not shown).

Sample12 16/34 16/34 Sample14 16/29 16/29 Sample16 16/36 16/36







Figure 6. Distribution of TOMM40 genotypes across 93 patient samples. Distribution of TOMM40 genotypes in a cohort of neuropathologically confirmed AD cases (AD cohort, n=64) and a control cohort without AD neuropathology (n=29) from the Center for Neurodegenerative Disease Research (CNDR) brain bank at the University of Pennsylvania. All were genotyped on the 3500xL instrument as described.

Conclusions

- The AmplideX[®] PCR/CE TOMM40 Kit (RUO) enabled single-base poly-T resolution, reproducible genotyping, and unambiguous data interpretation from cell-line, blood, buccal, and brain DNA.
- Matched specimens from blood and brain cells, and blood and buccal cells, demonstrated identical TOMM40 '523 genotypes.
- The combination of reproducible assay results and a simple, rapid, and scalable workflow may help standardize findings across many different sites in large-scale research studies.
- This assay has potential to improve risk assessments for LOAD, particularly when combined with APOE genotyping and/or other prognostic information.

- Crenshaw DG, et al. Using genetics to enable studies on the prevention of Alzheimer's disease. Clin Pharmacol Ther. 2013;93(2):177-85.
 Yu L, et al. *TOMM40* '523 variant and cognitive decline in older persons with *APOE E373* genotype. Neurology. 2017 Feb 14;88(7):661-668.



Equal contribution authors Research Use Only - Not for Use in Diagnostic Procedure Presented at AAIC 2017