A Unified PCR/CE Carrier Screening Workflow for CFTR, SMN1, and FMR1 Variants that Consolidates the Detection of SNVs, Indels, CNVs, and Triplet Repeat Expansions

Liangjing Chen, Elliot Hallmark, Brian C Haynes, Pranesh Rao, Shobha Gokul, Lawrence Manzano, Bradley Hall, and Gary J Latham Research and Development, Asuragen, Inc., Austin, TX

Summary

- Broad access to carrier screening requires streamlined, accurate, and scalable assays that can be run in different lab settings using readilyavailable instrumentation.
- We developed a PCR/CE assay* using AmplideX[®] chemistry that identifies 23 CFTR variants [1] recommended by ACMG/ACOG.
- The CFTR assay design is compatible with coinjection of SMN1 copy number and FMR1 CGG repeat-primed PCR amplicons.
- This assay approach enables a simple, singleplatform solution for the detection of multiple variant types across three of the most common carrier screening genes.







Introduction

About 5% of Western populations are genetic carriers for one of three disorders: Fragile X Syndrome (FXS), Spinal Muscular Atrophy (SMA), or Cystic Fibrosis (CF). Reliable detection of the pathogenic mutations that cause these disorders requires molecular assays that accommodate multiple classes of sequence variants, namely CGG repeat expansion in *FMR1* (FXS), copy number variants in *SMN1* (SMA), and single nucleotide variants (SNVs) and indels in *CFTR* (CF). These variants are each technically challenging to resolve, and current screening assays often rely on distinct technologies, platforms and workflows. We describe a simple and accessible multiplex PCR approach that can report genetic variants in *FMR1*, *SMN1*, and *CFTR* by capillary electrophoresis (CE) in a single workflow.

Materials and Methods

A prototype PCR/CE assay* was developed to detect and discriminate 23 mutations (SNVs and indels) in the *CFTR* gene. Multiplex primers were derived using a novel pipeline and designs that allowed detection of proximal variants. Primers were tagged with a NED dye and configured with different amplicon mobilities to optimize multiplex co-detection of *CFTR* mutations and resolve *SMN1* and *FMR1* PCR amplicons generated from AmplideX[®] PCR/CE *FMR1*^{**} and *SMN1*^{**} kits using a single CE injection on the Applied Biosystems[™] 3500 Series Genetic Analyzer.



Figure 2. Specific and Multiplexed Allelic Discrimination. Seven SNPs are located within <100 bp *CFTR* genomic region. Using a proprietary primer design strategy, specific differentiation of wt and mut was achieved across all seven mutations. **A)** Cell line NA12878 from the Coriell Cell Repository (CCR) was used for wt alleles (top). **B)** A gBlock bearing all seven mutant alleles shows specific detection of mutant alleles. **C)** A mixture of gBlock and cell-line DNA demonstrates multiplexed detection in a single PCR reaction.



Figure 3. Multiplex CFTR Assay Targeting 23 Mutations Resolved on the 3500xL CE. A) NA12878 manifests 23 wt peaks at the expected locations. B) Mixture of 23 gBlocks with mutant alleles for each of the 23 mutations displays peaks shifted 3 bp shorter than the corresponding wt peak size per the design. C) Cell line NA11280 with heterozygous mutations 621+1G>T and 711+1G>T. Both mutant alleles were observed in addition to the wt allele. D) Cell line NA11283 with heterozygous mutations A455E and dF508. Both mutant alleles were observed in addition to the wt allele.



Figure 4. Representative Electropherograms for Co-detected FMR1 CGG Repeats, SMN1 Copy Numbers and CFTR Mutations. All amplicons were pooled and co-injected with the ROX 1000 ladder from Asuragen. A) Cell line NA12878 B) Cell line NA11280 C) Cell line CD00014. To simplify comparison of electropherograms, the full-mutation pile-up peak for the 645 CGG trace is not shown. D) Cell line NA04025.



Figure 1. PCR/CE Primer Design Strategy Enables Differentiation and Detection of Mutant (mut) and Wild-type (wt) Alleles. The schematic shows the primer design strategy for PCR and a virtual electropherogram for an unspecified number of A) mutant and B) wildtype amplicon pairs. For the detection of multiple mutations (>50 in a single dye channel), the relative positions of the gene-specific (GS) common primer and allele-specific (AS) primers were adjusted to ensure complete separation between different variants on CE (panel C). Wt and mut peaks are separated by 3 nucleotides (spacer) on the electropherogram in our designs.

CCR	FMR1 CGG Repeat #		SMN1 Copies	CETP Mut Detected
	Allele 1	Allele 2	Smill Copies	CI IN MUL DELECLEU
NA12878	30		2	wt
CD0008	28	30	2	2184delA
NA00897	30		2	dF508
NA00999	30	56	2	621+1G>T
NA07464	29		2	R553X
NA07732	29		2	dF508
NA11277	32		2	dI507
NA11280	31	41	2	621+1G>T 711+1G>T
NA11283	30		2	A455E, dF508
NA11288	36		2	dF508
NA11496	29	33	2	G542X
NA11723	20	30	2	W1282X
NA11859	20	29	2	2789+5G>A
NA12585	30		2	R1162X
NA12960	30		2	R334W
NA13423	29		3	G85E
NA20737	29	38	2	R347H
NA21847	19	23	2	dF508
NA00232	55		0	wt
NA03815	20		1	wt
NA05756	30		2	wt
NA17117	21		3	wt
NA13664	30	52	2	wt
CD00014	56		2	wt
NA20242	30	74	2	wt
NA06891	119		2	wt
NA02039	20	>200	2	wt
NA04025	>200		2	wt

Table 1. Single-platform, Combined Detection of CFTR Mutations,

Conclusions

- We report a single-tube, PCR/CE assay that can resolve all 23 ACMG recommended *CFTR* variants.
- Our results demonstrate feasibility for a streamlined, single-platform, multi-modal, PCR-based assay that can unify the detection of *FMR1*, *SMN1*, and *CFTR* variants.
- This technology has the potential to address carrier screening in decentralized laboratory settings using standardized reagents, existing instrumentation, and fit-for-purpose analysis software.

Reference

 Watson MS, Cutting GR, Desnick RJ, Driscoll DA, Klinger K, Mennuti M, Palomaki GE, Popovich BW, Pratt VM, Rohlfs EM, Strom CM, Richards CS, Witt DR, Grody WW: Cystic fibrosis population carrier screening: 2004 revision of American College of Medical Genetics mutation panel. Genet Med 2004, 6:387e391.

SMN1 Copy Number, and FMR1 Triplet Repeats. Nineteen CFTR mutations (within the 23 panel) were detected in 17 mutant CFTR cell lines. All FMR1 CGG repeat numbers and SMN1 copy numbers were correctly called using AmplideX[®] PCR/CE Reporter software with FMR1 and SMN1 assay-specific modules. The co-detection results were also 100% concordant with separate CE injection for each individual assay.

*Product under development. Future availability and performance cannot be ensured. **Research Use only. Not for use in diagnostic procedures. Presented at ASHG 2018 -2892W

