Accurate Single-Tube Quantification of SMN1 and SMN2 **Copy Numbers Using a Rapid and Streamlined PCR/CE Assay Evaluated at Two Different Laboratories**

<u>Márcia E Oliveira</u>¹, Stela Filipovic-Sadic², Shobha Gokul², Justin Janovsky², Walairat Laosinchai-Wolf², John Milligan², Gary J Latham², and Rosário Santos¹

¹Molecular Genetics Unit, Centro Genética Médica Dr. Jacinto Magalhães, Centro Hospitalar Porto, Porto, Portugal; ²Asuragen, Inc., Austin, TX USA

Summary

- Spinal Muscular Atrophy (SMA) is an autosomal recessive neuromuscular disease and a leading genetic cause of infant mortality.
- The severity of SMA is often modulated by the copy number of the paralogous SMN2 gene.

Table 1. Phase 1 Results from Sites 1 and 2 for Normalized, Pre-binned SMN1 and SMN2 **Copy Numbers.** Both sites tested a set of 12 residual clinical samples spanning 0 to ≥ 4 copies for both SMN1 and SMN2 as part of Phase 1 training and proficiency. All results matched expected copy numbers determined using MLPA.

Sample Name	Expected SMN1 Copies	Site 1	Site 2	Expected SMN2 Copies	Site 1	Site 2
AS_CHP_I-01	0	0.0	0.0	2	1.9	1.9
AS_CHP_I-02	2	1.9	1.8	1	0.9	0.9
AS_CHP_I-03	2	2.1	1.8	1	1.0	0.9
AS_CHP_I-04	0	0.0	0.0	3	2.8	2.6
AS_CHP_I-05	0	0.0	0.0	3	2.8	2.8
AS_CHP_I-06	4	4.1	3.6	0	0.0	0.0
AS_CHP_I-07	3	2.9	2.6	1	0.9	0.9
AS_CHP_I-08	1	1.1	1.0	1	1.0	1.1
AS_CHP_I-09	1	1.0	0.8	2	1.9	1.8
AS_CHP_I-10	2	2.1	1.8	2	1.9	1.9
AS_CHP_I-11	0	0.0	0.0	4	3.8	3.5
AS_CHP_I-12	1	1.0	0.9	4	4.9	4.6

- We report a two-site evaluation of a single-tube SMN1/2 PCR/CE assay that quantifies SMN1 and SMN2 copy number (Site 1: Centro Hospitalar Porto, Site 2: Asuragen).
- Across both sites, the SMN1/2 PCR/CE assay accurately genotyped 60 samples, including 15 SMA and 13 carrier samples, using a rapid and simple workflow.

Introduction

Spinal Muscular Atrophy (SMA) is an autosomal recessive neuromuscular disease and a leading genetic cause of infant mortality. SMA is commonly caused by homozygous exon 7 deletions in the survival motor neuron 1 gene (SMN1). The severity of SMA is modulated by the copy number of the paralogous SMN2 gene. With significant progress made toward disease modifying treatments for SMA, such as the first approved drug nusinersen (SPINRAZA[®]), early detection of SMA along with knowledge of SMN2 copy number is critical. Here, we report the two-site evaluation of a single-tube assay that quantifies SMN1 and SMN2 copy number.

Materials and Methods

A prototype SMN1/2 single-tube PCR was developed from

Phase 1 and 2 Results



Figure 3. Site-to-Site Comparisons of SMN1 and SMN2 Normalized Copy Numbers from **Phase 1 and 2 Studies.** *SMN1* and *SMN2* plots of averaged, normalized copy numbers for 60 samples tested on the 3130 genetic analyzer. The site-to-site correlation coefficient was >0.98 for both SMN1 and SMN2 copy numbers.

AmplideX[®] PCR/CE SMN1 reagents* (Asuragen). Amplicons generated from whole blood genomic DNA were resolved by capillary electrophoresis (CE) on 3 different genetic analyzers: 3130, 3500, and SeqStudio (Thermo Fisher Scientific). Gene copy numbers were calculated from peak area ratios that were normalized to a plate calibrator and binned as 0, 1, 2, 3, or ≥ 4 copies. Two laboratories (one in Porto, Portugal and one in the US) evaluated 60 residual clinical samples, with reference genotypes obtained by MLPA (P021-A2, MRC-Holland). The two-site evaluation was executed in two phases with 12 samples used for training (Phase 1) and 48 for clinical evaluation (Phase 2) using the 3130 genetic analyzer. In addition, 60 samples, including a subset of 16 samples from the initial study, were used to compare the 3500 and SeqStudio CE instruments at Site 2. Data were analyzed using GeneMapper Software (Thermo Fisher Scientific) and copy numbers were calculated using AmplideX PCR/CE SMN1/2 Macro (Asuragen).

1A	PCR Cycles	CE Instrument, Capillary Length	Injection	Pre-Run	Run	Oven Temp
	27	3130xl, 36 cm	2.5 kV, 20 s	15 kV, 900 s	15 kV, 1500 s	60°C
	25	3500xL, 50 cm	2.5 kV, 20 s	15 kV, 900 s	19.5 kV, 2100 s	60°C
	25	SeqStudio	6 kV, 2 s	13 kV, 180 s	6 kV, 3000 s	60°C



Table 2. Site-to-Site Call Concordance and Comparisons with Reference Results. Binned copy-number calls for SMN1 were 100% concordant across both sites and with a Reference method. SMN2 calls were concordant for 59/60 (98.3%) samples between the two sites, for 60/60 (100%) samples between Site 1 and the Reference method, and for 59/60 (98.3%) samples between Site 2 and the Reference method.

		0 сору	1 сору	2 сору	3 сору	≥4 copy
SMN1	Site 1	15	13	14	12	6
	Site 2	15	13	14	12	6
	Reference	15	13	14	12	6
SMN2	Site 1	7	14	19	11	9
	Site 2	7	15	18	11	9
	Reference	7	14	19	11	9



Figure 4. CE Instrument Agreement for Normalized Copy Numbers. SMN1 and SMN2 plots of averaged, normalized copy numbers for 60 samples (16/60 samples from Phases 1 & 2 plus 44 additional blood derived samples tested at Site 2 only). 3500 and SeqStudio comparison demonstrated a correlation coefficient of >0.99 for both SMN1 and SMN2. Correlation coefficients were also >0.99 for all pairwise combinations of SMN1 and SMN2 copy-number calls using the 3130, 3500, and SeqStudio for 16 common samples.

0.300-0.600	0.305-0.615	1
0.700-1.125	0.700-1.070	2
1.170-1.570	1.095-1.650	3
≥1.615	≥1.675	≥4

Figure 1. PCR/CE Assay and Data Analysis Workflow. A) PCR cycles and injection conditions recommended for specific CE instruments. B) Data analysis workflow. Normalized ratios are calculated for each peak by dividing peak area by endogenous control peak area and normalizing to the same ratio derived from the calibrator sample (included in each batch). Normalized ratios are then placed within copy number bins, and results are reported as 0, 1, 2, 3, or \geq 4 copies.

Results



Figure 2. The SMN1/2 PCR/CE Assay Workflow is Streamlined from Sample-to-**Answer.** The assay can be performed with a turnaround time of ~3 hours for 1 capillary electrophoresis injection (4-24 samples) with ~50 min of total hands-on time.

*Research Use Only - Not for use in diagnostic procedures Presented at ESHG 2019

Conclusions

- The single-tube AmplideX[®] PCR/CE SMN1/2 assay* accurately quantifies 0, 1, 2, 3, or \geq 4 gene copies for both SMN1 and SMN2 from blood specimens using a ~3 hour workflow.
- Copy-number calls for SMN1 were 100% concordant across two laboratory sites and with a Reference method; SMN2 copy-number calls were 98% concordant across the sites (using a 3130 CE instrument).
- Accurate copy-number quantification was also observed on the 3500 and SeqStudio genetic analyzers, with instrumentto-instrument correlations of >99%.





