

A Multiplex PCR/CE *CFTR* Assay Resolves Zygosity of the 23 ACMG/ACOG-recommended *CFTR* Variants and Sizes poly-T and TG Repeats in a Single Tube

Pranesh Rao, Darshana Patel, Connor Parker, Gary J Latham, Liangjing Chen and Bradley Hall
Asuragen, Inc., Austin, TX

Summary

- Mutations in the *CFTR* gene on Chromosome 7 are responsible for Cystic Fibrosis (CF), which affects 1 in 3300 births.
- Reliable detection of CF mutations informs disease diagnosis, therapy decisions, carrier screening and pre-natal testing.
- We developed a PCR/CE assay that accurately genotypes all 23 ACMG/ACOG-recommended variants plus 7 additional high-prevalence variants and IVS8 (poly-T/TG) using a simple and rapid one-tube workflow.

Introduction

Cystic Fibrosis (CF) is an autosomal recessive condition caused by mutations in the CF transmembrane conductance regulator gene (*CFTR*) responsible for chlorine transport across the cell membrane. Mutations in the coding region result in dysfunctional ion transport¹. This typically results in increased chloride concentration in sweat, thicker mucus linings in bronchi, impaired pancreatic exocrine function and intestinal absorption.

ACMG and ACOG² recommends a panel of 23 mutations (CF23) on the basis of mutation frequency in CF patients in the US. Poly-T sizing, which affects splicing of exon 10 primarily for 5T alleles, is a modifier of R117H³ and recommended by ACMG as a reflex test. Further, expansion of a dinucleotide TG-repeat directly 5' of the poly-T tract in phase with a 5T repeat causes reduced active *CFTR* protein expression⁴ and should be accommodated within an assay design. Recently, Beauchamp et al.⁵ published NGS results of *CFTR* mutant allele frequency for an ethnically diverse US population of >115K individuals that includes mutations of higher frequency than many of those in the CF23 panel.

Here we describe a prototype assay for the recommended CF23 mutations plus 7 additional high-prevalence variant alleles that sizes/ phases the poly-T/TG tract in a single-tube PCR/CE assay.

Materials and Methods

Cell-line genomic DNA samples with known *CFTR* mutations were obtained from Coriell Cell Repository (CCR; n=62) or Asuragen (n=11). Genomic DNA was isolated from whole blood (WB) samples of 124 presumed healthy donors. Genomic DNA was amplified using AmpliDex[®] PCR chemistry, and resolved by capillary electrophoresis (CE) on a 3500xL genetic analyzer (Thermo Fisher Scientific) within eight hours from DNA to genotype. Alleles and zygosity were differentiated by mobility shift compared to a ROX ladder. Variants of the poly-T tract were further distinguished by color channel (NED, FAM or HEX). Sanger sequencing was used as a comparator method to confirm cell lines (n=73), WB carriers (n=4) and WB wild type (n=5).




Figure 1. Design of a Prototype Single-tube, Multiplexed *CFTR* PCR/CE Assay. A) Allele-specific primers amplify wild-type (WT) or mutant (MUT) allele variants with a common primer labeled with a NED fluorophore. A three-nucleotide (nt) spacer on the WT primer differentiates the allele by size. The poly-T IVS8 allele was sized by dye channel and the poly-TG track was sized by mobility shift. B) A representative CE trace for the 32-plex assay. Allele common names are listed below each peak. A heterozygous allele is observed by the presence of both the MUT and the WT peaks (i.e. R1162X). CF23 allele names are given in black, additional alleles are shown in orange.

Results




Figure 2. Zygosity is Accurately Resolved for Each Mutant *CFTR* Allele. CE traces were visualized by GeneMapper 5.0. Peaks were binned based on expected sizes. Amplicons from WT alleles were observed in teal bins and MUT in grey except for IVS8. A) NA12878 is homozygous WT for all alleles. This cell line does not have a 9T allele (in NED ch.). B) A mix of synthetic gBlock containing MUT alleles is MUT for all alleles. C) NA12585 is WT for all alleles except R1162X [ARG1162TER] where both WT and MUT alleles are observed, where the mutant allele peak is offset from its corresponding WT peak by 3bp. This would be common for a CF carrier. D) NA07441 is a compound heterozygote where both WT and MUT alleles are observed for 3120+1G>A and 621+1G>T variants. This cell line has a 9T allele phased with 10TG (yellow bin). E) NA01531 is a homozygous mutant for F508del while F) NA21847 is a carrier for F508del with a single affected allele. All samples with F508del were also observed to have at least one 9T allele for IVS8. Both D and E would be typical of a CF affected patient while C and F would reflect a CF carrier.




Figure 3. The Poly-T/TG Tract (IVS8) is Accurately Genotyped and Phased using a Novel Assay Design. Dye-specific peaks flag poly-T repeat length: 5T (HEX), 7T (FAM) and 9T (NED). Peak size (bp) determines the poly-TG size and phases both tracts. A) NA11723 has both 5T and 7T. The 5T phases in the 12TG bin and the 7T phases in the 10TG bin. B) NA11280 has both a 7T and 9T. The 7T is phased with 11TG and the 9T is phased with 10TG. The other peaks in the NED channel (A & B) correspond to the wt/mut peaks of the remaining 30 alleles C) Sanger sequencing was used to determine sizing truth for the IVS8 tract.

ASGN Sanger Sequencing	Concordance with comparator assay		ASGN 32-plex <i>CFTR</i> PCR/CE		Zygosity Concordance	
	homozygous	heterozygous	wt	mut	wt/mut	wt/mut
	wt	wt/mut	4288	-	-	4288/4288 (100%)
	mut	wt/mut	-	8	-	8/8 (100%)
	heterozygous	wt/mut	-	-	42/42	84/84 (100%)

Sample name	Poly-T	Poly-TG	Sample name	Poly-T	Poly-TG	Sample name	Poly-T	Poly-TG
NA11723	5T/7T	12TG/10TG	CD00003	7T/7T	11TG/11TG	NA06906	7T/7T	11TG/12TG
NA20242	5T/9T	11TG/10TG	CD00007	7T/7T	11TG/11TG	NA07464	7T/7T	11TG/12TG
NA13591	5T/9T	12TG/10TG	CD00014	7T/7T	11TG/11TG	NA12960	7T/7T	11TG/12TG
NA11277	7T/7T	10TG/10TG	NA12878	7T/7T	11TG/11TG	RU0005	7T/7T	11TG/12TG
NA11859	7T/7T	10TG/10TG	NA13033	7T/7T	11TG/11TG	RU0006	7T/7T	11TG/12TG
NA11860	7T/7T	10TG/10TG	NA13423	7T/7T	11TG/11TG	NA20837	7T/7T	11TG/12TG
NA12585	7T/7T	10TG/10TG	NA13654	7T/7T	11TG/11TG	CD00012	7T/7T	11TG/12TG
NA12926	7T/7T	10TG/10TG	NA20239	7T/7T	11TG/11TG	GMO4345	7T/7T	11TG/12TG
NA07857	7T/7T	10TG/10TG	NA20741	7T/7T	11TG/11TG	NA09999	7T/7T	11TG/12TG
NA12444	7T/7T	10TG/10TG	NA20836	7T/7T	11TG/11TG	NA11278	7T/7T	11TG/12TG
NA12785	7T/7T	10TG/10TG	NA20915	7T/7T	11TG/11TG	NA11284	7T/7T	11TG/12TG
NA12961	7T/7T	10TG/10TG	NA20925	7T/7T	11TG/11TG	NA11370	7T/7T	11TG/12TG
NA18801	7T/7T	10TG/10TG	RU0002	7T/7T	11TG/11TG	NA18803	7T/7T	11TG/12TG
NA18802	7T/7T	10TG/10TG	RU0007	7T/7T	11TG/11TG	NA21551	7T/7T	11TG/12TG
NA18803	7T/7T	10TG/10TG	RU0009	7T/7T	11TG/11TG	NA21847	7T/7T	11TG/12TG
NA20737	7T/7T	10TG/10TG	RU0010	7T/7T	11TG/11TG	NA04025	7T/7T	11TG/12TG
NA21069	7T/7T	10TG/10TG	RU0011	7T/7T	11TG/11TG	NA08338	7T/7T	11TG/12TG
RU0003	7T/7T	10TG/10TG	NA06891	7T/7T	11TG/11TG	CD00008	7T/7T	11TG/12TG
RU0008	7T/7T	10TG/10TG	RU0004	7T/7T	11TG/11TG	CD00009	7T/7T	11TG/12TG
						NA08897	7T/7T	11TG/12TG
						NA07441	7T/7T	11TG/12TG
						NA07732	7T/7T	11TG/12TG
						NA11275	7T/7T	11TG/12TG
						NA11280	7T/7T	11TG/12TG
						NA11288	7T/7T	11TG/12TG
						NA11472	7T/7T	11TG/12TG
						NA18800	7T/7T	11TG/12TG
						NA2924	7T/7T	11TG/12TG
						NA21080	7T/7T	11TG/12TG
						NA22063	7T/7T	11TG/12TG
						RU0001	7T/7T	11TG/12TG
						NA01531	9T/9T	10TG/10TG
						NA11283	9T/9T	10TG/10TG
						NA11496	9T/9T	10TG/10TG
						NA18886	9T/9T	10TG/10TG

Figure 4. The Prototype Assay was 100% Concordant with Comparator Method. Cell lines (n=73) were sequenced to determine truth for each of the 32 alleles covered in the assay. A) The assay was 100% concordant for both WT and MUT alleles observed (26 mut cell lines + 6 mut gBlocks) and also to published data. B) Poly-T/TG sizing and phasing was 100% concordant with comparator for each of the cell lines and with published results for IVS8 where available (n=43)⁶.

#RU000x: Asuragen's internal cell lines.
[†G85C, N1303X, P5L, S945L, R1070Q and L206W]: No Coriell cell line DNA available for these variants.




Figure 5. Analysis of 124 Presumed Normal WB Donors Identified Four Mutant Alleles. Of the four carriers, three had an F508del mutant allele (A) which is the most common *CFTR* mutation. B) One carrier had an R117H mutant allele. All results were confirmed by sequencing. A carrier frequency of one in 31 (3.2%) was observed for the 32 allele set which is similar to published carrier rate observed in the US population (~3% or one in 35)⁷. C) The genotypic distribution of poly-T/TG suggests 7T/7T and 11TG/11TG are most common. All the F508del carriers (3/3) also had a 9T allele which corroborates the observations from literature³. The donor with a heterozygous R117H mutant allele did not have a 5T.




Figure 6. The Prototype 32-plex Assay Covers 86% of Pathogenic Alleles as Reported by Beauchamp et al. A) Allele frequency by Beauchamp et al. (2019) or CFTR2 covered by various on-market *CFTR* kits compared with the prototype ("Asuragen Prototype"). This prototype is in active development towards 61 variants ("Asuragen Dev") that would cover >92% allele frequency based on a more ethnically diverse population reported by Beauchamp et al. B) compared to CFTR2⁸ C).

Conclusions

- We developed and evaluated a 32-variant PCR/CE *CFTR* assay that called 4672 *CFTR* (wt+mut) alleles from 73 cell lines with 100% agreement to reference genotypes.
- The assay accurately resolved SNVs, indels and repetitive sequences, and included a novel strategy for dye-specific coding and fragment sizing to genotype poly-T/TG.
- Genotyping was accomplished from a single input, PCR and CE injection using a workflow that readily slots within an 8-hour shift, and with the flexibility to support high-throughput screening or lower-volume diagnostic applications.
- The current panel content covers 86% of the pathogenic variants reported in the US population and all mutations recommended by ACMG/ACOG guidelines. This coverage is similar to existing *CFTR* assay kits, but has a simpler procedure, fewer handling steps, and a single-tube configuration.

*This product is under development. Future availability and performance cannot be ensured. Presented at AMP 2019-G022

