


# Analytical Validation of a Multiplex PCR/CE Assay for Simultaneous Determination of *SMN1*/*SMN2* Exon 7 Copy Number and Variant Status

## Summary

- Spinal Muscular Atrophy (SMA) is an autosomal recessive neuromuscular disease that results from mutation of the survival motor neuron 1 gene (*SMN1*), where disease severity is modulated by the *SMN2* copy number.
- The AmpliDeX® SMA Plus Kit\* is a one-tube PCR/CE assay that quantifies *SMN1* and *SMN2* copy numbers and genotypes gene duplication “silent carrier” markers (c.\*3+80T>G and c.\*211\_\*212del) and a disease modifier variant (*SMN2* c.859G>C).
- Includes reagents, controls, and automated, push-button data reporting via AmpliDeX Reporter software that reports results in under 5 minutes, with ~60 min of total hands-on time (<4 hours sample to result).
- We validated the kit by testing 122 unique samples to assess single-site precision, accuracy compared to reference methods, analytical specificity, DNA input range, reagent stability, and other factors



AmpliDeX SMA Plus Kit\* Example Electropherogram Readout. Gene-specific peak area ratios relative to EC are normalized to a calibrator to produce normalized ratios, which are binned to determine copy number.

# AmplideX<sup>®</sup> SMA Plus Kit\*

Complete testing solution offering comprehensive analysis of *SMN1* and *SMN2*

## AmplideX<sup>®</sup> SMA Plus Kit

### Kit Components

2X PCR Mix	SMN Control
SMN1/2 HEX Primer Mix	Diluent
SMN Calibrator	ROX 1000 Size Ladder



**REACTIONS PER KIT**  
50, 100



**SAMPLE TYPES**  
Blood



**DNA INPUT**  
20-70 ng SeqStudio/3500/3730 Genetic Analyzers  
20-60 ng 3130 Genetic Analyzers



**CE INSTRUMENTATION**  
Applied Biosystems<sup>®</sup> Genetic Analyzers



**TIME-TO-RESULT**  
<4 HRS



**HANDS-ON-TIME**  
60 MINS

# AmplideX<sup>®</sup> PCR/CE Reporter Software

Simple and streamlined sample reporting at the push of a button



## FEATURES INCLUDE

- Quickly determine *SMN1* and *SMN2* copy number and detect silent carrier and disease modifier variants
- Batch-level QC and run metrics
- Project Management tool
- Export results in PNG and LIMS-compatible formats (Excel/CSV)

ILLUSTRATIVE

# Comparison to Multiple Common Reference Methods Shows Consistent Copy Number Results for *SMN1* and *SMN2*

		Consensus <i>SMN1</i> Copy Number (MLPA + ddPCR)					
		0	1	2	3	4	Sum
Measured <i>SMN1</i> Copy Number (AmplideX)	0	15	0	0	0	0	15
	1	0	16	0	0	0	16
	2	0	0	60	1	0	61
	3	0	0	0	20	0	20
	4	0	0	0	0	10	10
	Sum	15	16	60	21	10	122

***SMN1*: 99.2% (121/122)**

		Consensus <i>SMN2</i> Copy Number (MLPA + ddPCR)					
		0	1	2	3	4	Sum
Measured <i>SMN2</i> Copy Number (AmplideX)	0	25	0	0	0	0	25
	1	0	40	0	0	0	40
	2	0	0	33	0	0	33
	3	0	0	1	15	0	16
	4	0	0	0	1	6	7
	Sum	25	40	34	16	6	121

***SMN2*: 98.3% (119/121)**

Target	Percent Agreement
c.*3+80T>G	100% (132/132)
c.*211_*212del	100% (118/118)
c.859G>C	100% (132/132)

- 134 DNA samples (116 unique specimens, 7 unique isolation methods) measured in singleton on a 3500xL Genetic Analyzer and analyzed with AmplideX Reporter Software
  - Precision QC failures (PR, n=4 for *SMN1* and *SMN2*) were excluded from analysis
- Two reference methods were used to measure *SMN1* and *SMN2* exon 7 copy numbers (MLPA and ddPCR)
  - Samples with disagreement or ambiguous results in either method were excluded from analysis
- **Results were over 98% concordant to reference results for both *SMN1* and *SMN2***
  - Fewer QC failures observed in AmplideX® SMA Plus Kit\* than in ddPCR, which caused majority of ambiguous reference method results (n=6 of 9 for *SMN1*, n=8 of 10 for *SMN2*)

# Single-site Precision Shows Low Variability Across Multiple Runs, Days, Reagent Lots, and Operators

Sample	Sample Type	SMN1 Copies	SMN1 Copy Number Percent Agreement	SMN1 %CV (NR)	SMN2 Copies	SMN2 Copy Number Percent Agreement	SMN2 %CV (NR)
1	Blood	0	100% (40/40)	NA	2	100% (40/40)	8.0
2	Blood	1*	97.5% (39/40)	5.8*	3	96.8% (30/31)	5.8
3	Blood	1	100% (40/40)	6.3	4	100% (38/38)	5.8
4	Blood	4	100% (40/40)	4.5	1	100% (40/40)	6.3
5	Cell Line	0	100% (40/40)	NA	2	100% (40/40)	5.4
6	Cell Line	1	100% (40/40)	6.7	3	97.4% (37/38)	6.9
7	Cell Line	3	100% (38/38)	5.6	0	100% (39/39)	NA
8	Cell Line	2	100% (40/40)	7.9	2	100% (40/40)	8.4
9	Cell Line	2	100% (40/40)	6.9	2	100% (40/40)	5.8
		<b>Total:</b>	<b>99.7% (357/358)</b>		<b>Total:</b>	<b>99.4% (344/346)</b>	

Target	Percent Agreement
c.*3+80T>G	100% (279/279)
c.*211_*212del	100% (279/279)
c.859G>C	100% (279/279)

\*copy number and NR measurements based on SMN1 hybrid

- Study consisted of 9 DNA samples tested in duplicate in 20 batch runs over 10 days using two operators, three reagent lots, two thermal cyclers on a 3500xL Genetic Analyzer
  - QC failures were excluded from analysis
- Percent CV for normalized ratios (NR) and percent agreement with the mode were calculated by sample for both SMN1 and SMN2
- **Across all measurements, copy numbers were over 99% concordant for both SMN1 and SMN2**

# Results Are Consistent Across Full DNA Input Range Between Reagent Lots and CE instruments for All Relevant Genotypes

DNA Input	SMN1 Copy Number Percent Agreement	SMN2 Copy Number Percent Agreement
10	96.8% (121/125)	100% (123/123)
20	100% (128/128)	100% (128/128)
50	99.2% (126/127)	98.4% (124/126)
60	99.2% (123/124)	98.4% (124/126)
70	100% (124/124)	98.3% (118/120)
80	100% (89/89)	100% (95/95)
<b>Total</b>	<b>99.2% (711/717)</b>	<b>99.2% (712/718)</b>

Target	Percent Agreement
c.*3+80T>G	100% (735/735)
c.*211_*212del	100% (735/735)
c.859G>C	100% (735/735)

- Study consisted of 8 DNA samples tested in duplicate at 10, 20, 50, 60, 70, and 80 ng total DNA input using two reagent lots on the following Genetic Analyzers: 3500xL, 3730xl, 3130xl, and SeqStudio
- Sample panel covered 0, 1, 2, 3 copies of *SMN1*, 0, 1, 2, 3, ≥4 copies of *SMN2*
- Assay DNA Input range: 20-70 ng (3500xL, 3730xl, SeqStudio), 20-60 ng (3130xl)
  - QC failures (all instruments) and 80 ng measurements (3130xl only; above intended input range) were excluded from analysis
- Percent agreement with expected copy number results were calculated for each DNA input for both *SMN1* and *SMN2*
- **Across all measurements, copy numbers were over 99% concordant for both *SMN1* and *SMN2***

# Analytical Validation of a Multiplex PCR/CE Assay for Simultaneous Determination of *SMN1*/*SMN2* Exon 7 Copy Number and Variant Status

## Overall Performance and Conclusions

CE Platform	<i>SMN1</i> Copy Number Percent Agreement	<i>SMN2</i> Copy Number Percent Agreement	c.*3+80T>G Percent Agreement	c.*211_*212del Percent Agreement	c.859G>C Percent Agreement
3500xL	99.5% (739/743)	99.2% (723/729)	100% (781/781)	100% (767/767)	99.9% (780/781)
3730xl	100% (165/165)	100% (161/161)	100% (191/191)	100% (191/191)	100% (191/191)
3130xl	97.5% (154/158)	99.3% (152/153)	100% (191/191)	100% (191/191)	100% (191/191)
SeqStudio	100% (167/167)	100% (168/168)	100% (192/192)	100% (192/192)	100% (192/192)
<b>Total</b>	<b>99.4% (1225/1233)</b>	<b>99.4% (1204/1211)</b>	<b>100% (1355/1355)</b>	<b>100% (1341/1341)</b>	<b>99.9% (1354/1355)</b>

- Across all analytical validation studies, percent agreement was over 99% for *SMN1* and *SMN2* copy numbers and all variants detected
  - Including studies shown here as well as reagent stability and thermal cycler equivalency studies not shown
  - Consensus between ddPCR and MLPA as reference values for *SMN1* and *SMN2*, Sanger sequencing as reference for variants
- Analytical specificity (exclusivity) from method comparison, precision, and DNA input studies for samples with 0 copies of *SMN1* or *SMN2* and for all variants detected by the assay was 100%
- The expanded content of the AmpliX SMA Plus Kit\* combined with its convenience and versatility may provide an all-in-one option for laboratories interested in SMA carrier screening and diagnostics.
- Visit <https://asuragen.com/eshg-20/> to download a copy of this presentation