A Streamlined PCR/Nanopore Sequencing Carrier Screening Panel for Cystic Fibrosis, Spinal Muscular Atrophy, and Fragile X Syndrome

Jon Kemppainen, Pranesh Rao, Julie R Thibert, Walairat Laosinchai-Wolf, Ryan Routsong, Ted Markulin, Monica Roberts, Bradley Martin, Bryan Killinger, Chris Fraher, Gary J Latham, and Bradley Hall Asuragen, a Bio-Techne Brand, Austin, TX

Summary

- Cystic Fibrosis (CF), Spinal Muscular Atrophy (SMA), and Fragile X Syndrome (FXS) are three of the most common inherited genetic disorders, each with high carrier rates that often require distinct genotyping methods.
- We developed a prototype assay comprised of novel PCR enrichment and Nanopore sequencing to simultaneously detect SNV/indels, copy number variation, and repeat sizing in a unified workflow, without the need for manual analysis.
- The assay utilizes amplicon read depth and machine learning models to automate and streamline identification of key genetic variants specific to each disease.
- Assay performance was evaluated with a cohort of 94 cell-line and 207 whole blood samples for all three genes resulting in 96-100% accuracy for each variant class.

Introduction

Cystic Fibrosis (CF), Spinal Muscular Atrophy (SMA), and Fragile X Syndrome (FXS) are three of the most common inherited genetic disorders, each with high carter rates (-1/30 for CF, -1/50 for SMA and -1/250 for FXS). Historically, screening has required separate assays for each gene, often with Inefficient workflows and uneven detection rates across ethnicities. Additionally, genes for SMA and FXS (SMN)⁷ and FMR1, respectively) are technically challenging to characterize due to homologuous or GC-rich sequences. Repetitive oplymorphisms, exon detelbions and other CFTR variants can also be problematic to resolve.

Here, we demonstrate a prototype long-read sequencing (LRS) trio assay on Oxford Nanopore Technologies' (ONT) MinION platform. The assay identifies pathogenic variants across all exons of CFTR, exons 3-8 in SMN1/2, and sizes the CGG repeats in FMR1 with phased AGG status.

Materials and Methods

Cell-line genomic DNA samples were obtained from Coriell Cell Repository including HapMap controls, 0 to 24 copies of SM/H12, 48 unique pathogenic SNVsindels in CFTR representing >87% carrier prevalence' in the US population, and all major repeat expansion categories for FMRT, Cenomic DNA was isolated from whole blood of presumed healthy donors. DNA was amplified in two PCR enrichment reactions, barcoded, pooled, preped by ligation sequencing kt (ONT), and run on R9.4.1 flow cells using the MkTB (ONT). Data were analyzed using custom-developed software. The AmplideX* PCRV CE CFTR ktr, PCRVCE SM/H2 Plus ktr, PCRVCE FMRT ktr, Xpansion Interpreter-, Sanger sequencing, and 1000 Genomes were used to provide reference genotypes where appropriate.



Figure 1. PCRNRanopore Carrier Screening Panel Design and Workflow Identifies Pathogenic Variants for Three of the Most Common Inherited Genetic Diacoders in a Single Workflow. The combination of AmplideX PCR lachnology and Nanopore sequencing enables detection of highly GCrich repeat sequences in *KMR*; detection of AGO interrupts in *KMR*; highly homologous sequences such as *SMN* or *SMN2*, and pathogenic variants (SNVs and indeb) across all exons of *CFTR* and *SMN12*. One cellline sample was included for *FMR*? To yas an edge case for 7900 CGGs. Whole blood samples (re20) were included for *FMR*? and *SMN12* only since *CFTR* orthogonal data was unknown. Highlights of the streamlined workflow are shown under the paphic.



Figure 2. Sequence Data Differentiates SMM1 and SMM2 and Informs Silient Carrier Status and Disease Severity. A) Differentiation and assignment of reads to SMM1 and SMM2 by three different parado-specific variants align reads to each gene and inform copy number variation prediction model. B-D) Silent carrier (SC) or disease modifier (OM) variant alignments.

Cell Line





Figure 3. SMN12 PCRNanopore Assay Classifies Carrier, Silent Carrier, and Dieasee Modifier Status with 96-100% Accuracy. Colling accuracy for SMN1 and SMN2 corp numbers in a) 96 coll-line and 8) 227 whole blood samples. Hyperparameters for the decision tree model were selected using an 2020 traintist sell in a statified anothy selected for hold corps wildlick in scheme on an independent of 2020 traintist sell in a statified anothy selected for hold corps wildlick in scheme on an independent (crarge dashed culline). (2) Variant allele frequency at each SNV associated with silent carrier variant SSC (cr3+807-5, Hp nene), silent carrier variant SSC (cr2+12122de, independent) and disease modifier mutation DM (c.859G-5c, nght panel) were used to predict status of silent carrier or disease modifier acarch using F1 accose determined the optimal allele frequency threshold, with samples above the threshold being disalified as "Predicted Positiver" and samples below the threshold as "Predicted Negativer". The applichtm accuracity lottified SC1 (cr2) and DM mutation (cr3) "thresholding evidence of the encodes and 316231 (86.4%), SMM2 sample corp number agreed with the orthogonal method or corp number (CN) between 0 and >-3 gane copies.



Figure 4. Sequence Data Reveals Diverse Pathogenic Variants Across CF7R Exons. A) Heterocrypous FS08del, a 3-by deletion at GRCNa5 Acr7117.555.95 (a settion of phreyliamine at codon 508. B) Homocrypous missense variant M1101K, a T-A transversion at position chr7117.557.806 resulting in a protein change from methionen to sylane. C) Heterocrypous missense variant GS510, a GA-transition at position chr7.17.557.806 resulting in a protein change from glycine to aspartiz acid. D) Heterocrypous fiscation chr7.17.577.807.1778221 (Section 1778221) Feature in a female from a few by deletion GAA and insertion of G.

Table 1. CFTR Sample Agreement with Orthogonal Data for 94 Cell-line Samples (top) and 207 Presumen Normal Whole Blood Samples (bottom) Using Claid" or Coverage Differences for Exon Del/Dup Variant Identification. Assay detected 47/48 (07.9%) unique variants including two del/ou (CFTRediel 2.2, The Samples Called Samples Called Samples Called works and the Samples Called Samples Called Samples Called Works and Works and CFTRediel 9.2) unique variants including two del/ou samples Called Called Called Samples Called Samples Called Works and Works and CFTRediel 9.2) unique variants including two del/ou provents (20850) (19.0%) samples cared with orthogonal data.

Cell Line						
Sample level accuracy	Allele1/ Allele2	WT/WT	MUT/WT	MUT/MUT	Genotype agreement	
CFTR PCR/ Nanopore	WT/WT	48	2	0	48/50 (96.0%)	
	MUT/WT	0	20	0	20/20 (100%)	
	MUT/MUT	0	0	24	24/24 (100%)	

Whole Blood			nplideX PCR/CE CF				
Sample level accuracy	Allele1/ Allele2	WT/WT	WT/MUT	MUT/MUT	Genotype agreement		
CFTR PCR/ Nanopore	WT/WT	200	0	0	200/200 (100%)		
	WT/MUT	04	7	0	7/7 (100%)		
	митлиит	0	0	0	0/0 (100%)		

Table 2. FMRT categorical Agreement with Orthogonal Data for 95 Cell-Line and 227 Whole Blood Samples. ACMC categorical genotype boundaries are included for reference. All samples for White expected categories based AmplideA PCRUCE FMRT precision metrics (± 1: 0-70 repeats, ± 3: 71-120). Expanded samples to 16 40 CG2 (± 0.05%) amples and 454462 (± 0.5%) index. In // Orthogonal, CGG samples, the algorithm whole precision for 31 2222 (± 0.05%) samples and 454462 (± 0.5%) index. In // Orthogonal, CGG samples, the algorithm the algorithm definition a any videa videntified and more categories and earlies.

Sample Level gorical Accuracy	Normal <45 CGG	Intermediate 45-54 CGG	Premutation 55-200 CGG	Mutation >200 CGG		
Training	61	18	55	19	100%	100%
Cell line	72	5	14	4	100%	100%
Whole Blood	220	7	0	0	100%	100%



Figure 5. Predicted Risk of FMR1 Expansion Based on AGG Interruption Status. A control of 26 intermediate and premutation alleles were assessed using Asuragen Xpansion Interpreter (X) and PCR/ Nanopore using a custom algorithm. Genotypes were in 100% agreement with X1 for the absolute number of AGG interruption within the CGG repeat.

Table 3. Seventeen Carriers (8.2%), Including Those with Intermediate or Larger FMRT Expansions, were Identified in a Presumed Normal Cohort (n=207) Using the PCRNanopore Panel. One donor sample (SID530 in bold) was a carrier for both CFTR and SMN1. All variants were confirmed by orthogonal assay.

Sample ID		CFTR Allele 1	CFTR Allele 2						
SID537	46	WT	WT	2	0	(-)	(-)	(.)	FMR1 Intermediate
SID542	47	WT	WT	2	0	(-)	(-)	(.)	FMR1 Intermediate
SID497	48	WT	WT	2	1	(-)	(-)	(-)	FMR1 Intermediate
SID428	30,54	WT	WT	2	2	(-)	(-)	(.)	FMR1 Intermediate
SID535	39,47	WT	WT	2	2	(-)	(-)	(.)	FMR1 Intermediate
SID611	24,46	WT	WT	3 EXP 2 PRED	2	(-)	(-)	(-)	FMR1 Intermediate
SID481	31	F508del	WT	2	1	(-)	(-)	(.)	CFTR
SID562	30	F508del	WT	2	2	(-)	(-)	(-)	CFTR
SID569	30	F508del	WT	2	2	(-)	(-)	(.)	CFTR
SID637	30	F508del	WT	2	3	(-)	(-)	(-)	CFTR
SID489	30	R117H	WT	2	1	(-)	(-)	(.)	CFTR
SID420	31	R74W (VUS)	WT	3	1	(-)	Positive	Positive	CFTR
SID530	30	F508del	WT	2	1	(-)	Positive	Positive	CFTR, SMN1 SC
SID461	20,30	WT	WT	1	1	(-)	(-)	(-)	SMN1
SID546	29	WT	WT	1	2	(-)	(-)	(-)	SMN1
SID578	20,33	WT	WT	2	2	(-)	Positive	Positive	SMN1 SC
SID589	30	WT	WT	2	2	(-)	Positive	Positive	SMN1 SC

Conclusions

- The prototype PCR/Nanopore assay accurately resolves multiple challenging variants for three of the most common carrier screening genes using a streamlined workflow.
- At least 96 samples can be combined in a single run using only two PCR enrichment reactions. This has potential to reduce carrier screening turn around times when paired with other similar assays (see Poster # Q055).
- FMR1 CGG repeats and AGG interruptions associated with expansion risk from mother-to-child were detected in the same workflow.
- Detection of a dual carrier (CFTR and SMN1) highlights the importance of a unified carrier screening approach.
- In over 300 samples tested across each gene, the PCR/Nanopore assay agreed with the orthogonal method for SMN1 SNVs/indels (>99% of samples), SMN1 CN (96.6%), SMN2 CN (98.4%), CFTR SNVs/indels/ exon deletions (99.0%), FMR1 repeat categories (100%) and FMR1 AGG number (100%).

References

 Beauchamp, K. A. et al. Sequencing as a first-line methodology for cystic fibrosis carrier screening. Genet Med 37, 773–8 (2019).
Zheng, Z.: et al. Symphonizing pileup and full-alignment for deep learning-based long-read variant calling. Bioxiv 2021.12.29.474431 (2021).



Poster # G057

"This product is under development. Future availability and performance to be determined. For Research Use Only. Not for Use in Diagnostic Procedures. "Xpansion Interpreter" is a laboratory-developed test. Analytical and clinical performance have not been reviewed by the FDA All authors Thate the financial relationship to disclose: Employment by Asuragen.