An Amplification-based, Nanopore Carrier Screening Panel Resolves Clinically-Relevant Variants in CFTR, SMN1/2, HBA1/2, HBB, and FMR1 in a Unified Workflow

Bradley Hall, Bryan Killinger, Christopher J. Fraher, Bradley Martin, Monica Roberts, Jonathan Turner, Pranesh Rao, Ryan Routsong, Gary J. Latham Asuragen, a Bio-Techne Brand, Austin TX

G R

DM; c.859G>C

c

SMN2

Exon8

-2-

SUMMARY

- Cystic Fibrosis (CF), Spinal Muscular Atrophy (SMA), Hemoglobinopathies, and Fragile X Syndrome (FXS) are among the most commonly inherited genetic disorders, each with high carrier rates that often require distinct genotyping methods.
- We explored whether novel PCR enrichment and nanopore sequencing could detect multiple classes of variants across these genes in a single workflow.




Table 1. Pathogenic and likely pathogenic HBA1/2 and HBB variants detected using the PCR/nanopore assay. 145 CL and 247 WB samples were tested with the combined HBA1/2 and HBB assay. Of these, 102 CL and 235 WB samples were genotyped as wild-type (aa/aa). Sample variants are heterozygous (HET) unless listed as homozygous (HOM) or compound heterozygous (COMP). COMP included SEA/anti3.7 (HG01281), FIL/SEA (NA10797). Positive samples are all cell line unless indicated.

Gene	Variant	Туре	Positive Samples	Called Correctly 100%	
HBA	SEA	α ⁰ -thal HBA1/2 del	11 HET 1 HOM 2 COMP		


Table 3. FMR1 Categorical Agreement with Orthogonal Genotypes for 99 Cell-Line and 232 Whole Blood Samples. ACMG categorical genotype boundaries are included for reference. All samples fell within expected categories based AmplideX PCR/CE FMR1 precision metrics (± 1: 0-70 repeats, ± 3: 71-119). All expanded samples, including full mutations up to 940 CGG repeats, were flagged correctly. Additionally, CGG sizing was accurate within precision for 321/331 (97.0%) samples and 443/454 (97.5%) alleles. In 7/11 samples, the algorithm accurately called one of two alleles when two similar sized alleles (1-3 repeats difference) were present. The algorithm identified a previously unidentified minor mosaic allele in the remaining three samples that did not change the categorical call

- The assay utilizes PCR enrichment, coverage-based copy number, and machine learning models to automate and streamline identification of SNVs, INDELs, Exon del/dups, SVs, CNVs and STRs specific to each disease.
- The assay was optimized with 168 cell line samples and independently evaluated with 249 whole blood samples across the seven genes to identify potential carriers from presumed normal donors.


INTRODUCTION


Cystic Fibrosis (CF), Spinal Muscular Atrophy (SMA), Fragile X Syndrome (FXS), and Hemoglobinopathies, including alphaand beta-thalassemias, are each within the 10 most prevalent inherited diseases by carrier rate. Pathogenic variants in the associated genes can be difficult or intractable to identify using short-read sequencing alone. As a result, genotyping often relies on multiple non-NGS methods, highly trained (2A) personnel, and sophisticated or manual data analysis to resolve. Even still, uneven detection rates across ancestries persist.

Figure 1. PCR/nanopore Carrier Screening Panel Design and Workflow Identifies Pathogenic Variants for 5 of the Most Common Inherited Genetic Disorders in a Single Workflow. The combination of AmplideX® PCR technology and Nanopore sequencing enables detection of short tandem repeats (STRs), including the highly GCrich in FMR1, PolyT/TG region of CFTR, and detection of AGG interruptions in FMR1; differentiates and quantitates highly homologous sequences such as SMN1/2 or HBA1/2; and identifies pathogenic or informative variants (SNVs and indels) across CFTR, HBA1, HBA2, HBB, and SMN1/2. The majority of samples were tested across the full panel, especially for whole blood. For a subset of samples, only specific gene data was analyzed and compared, especially if a variant was known in the gene of interest. Highlights of the streamlined workflow are shown under the graphic.



HBA 3.7del α ⁺ -thal HBA2 del 14 HET CL, 4 HET WB 3 HOM CL 100% HBA anti-3.7 HBA2 del 5 HET CL, 8 HET WB 1 COMP CL 100% HBA 4.2del HBA2 del 2 100% HBA 4.2del HBA2 del 2 100% HBA1 c.427T>C (Hb Constant Spring) SNV/indel 1 1 HBA1 c.207C>G (Hb G-Philadelphia) SNV/indel 1 1 HBB (Hb Lepore-Baltimore) Exon 1, 2 del 2 100% HBB C.19G>A (Hb C) SNV/indel 1 100% HBB c.19G>A (Hb C) SNV/indel 2 100% HBB c.126_129del SNV/indel 3 100% HBB c.126_129del SNV/indel 2 100% HBB c.137C>G SNV/indel 1 100% HBB c.138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100%	HBA	FIL	αº-thal HBA1/2 del	1 HET 1 COMP	100%
HBA 3.7 del HBA2 del 3 HOM CL 100% HBA anti-3.7 HBA2 dup 5 HET CL, 8 HET WB 1 COMP CL 100% HBA 4.2del HBA2 del 2 100% HBA 4.2del HBA2 del 2 100% HBA 4.2del HBA2 del 2 100% HBA2 c.427T>C SNV/indel 1 1 HBA1 C.427T>C SNV/indel 1 1 HBA1 C.207C>G SNV/indel 1 1 HBB4 (Hb Censtant Spring) Exon 1, 2 del 2 100% HBB (Hb Lepore-Baltimore) Exon 1, 2 del 2 100% HBB c.19G>A (Hb C) SNV/indel 1 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.126_129del SNV/indel 3 100% HBB c.126_129del SNV/indel 2 100% HBB c.137C>G SNV/indel 1 </th <th>HBA</th> <th>THAI</th> <th></th> <th>1</th> <th>100%</th>	HBA	THAI		1	100%
HBAanti-3.7HBA2 dup1 COMP CL100%HBA4.2delHBA2 del2100%HBA2c.427T>C (Hb Constant Spring) (Hb G-Philadelphia)SNV/indel1HBA1c.207C>G (Hb G-Philadelphia)SNV/indel1HBB(Hb Lepore-Baltimore)Exon 1, 2 del2100%HBB(Sicilian [deltabeta]0-Thal)Exon 1-3 del1100%HBBc.19G>A (Hb C)SNV/indel2100%HBBc.20A>T (Hb E)SNV/indel5100%HBBc.20A>T (Hb S)SNV/indel3100%HBBc.126_129delSNV/indel2100%HBBc.93-21G>ASNV/indel1100%HBBc137C>GSNV/indel1100%HBBc.316-197C>TSNV/indel1100%HBBc.82G>T (Hb Knossos)SNV/indel1100%	HBA	3.7del			100%
c.427T>CSNV/indel1HBA2(Hb Constant Spring) c.207C>GSNV/indel1HBA1(Hb G-Philadelphia)SNV/indel1HBB(Hb Lepore-Baltimore)Exon 1, 2 del2100%HBB(Sicilian [deltabeta]0-Thal)Exon 1-3 del1100%HBBc.19G>A (Hb C)SNV/indel2100%HBBc.20A>T (Hb S)SNV/indel2100%HBBc.20A>T (Hb S)SNV/indel3100%HBBc.126_129delSNV/indel2100%HBBc.93-21G>ASNV/indel1100%HBBc137C>GSNV/indel1100%HBBc.316-197C>TSNV/indel1100%HBBc.82G>T (Hb Knossos)SNV/indel1100%	HBA	anti-3.7	HBA2 dup		100%
IBA2 (Hb Constant Spring) c.207C>GSNV/indel1IBA1(Hb G-Philadelphia)1IBB(Hb Lepore-Baltimore)Exon 1, 2 del2IBB(Sicilian [deltabeta]0-Thal)Exon 1-3 del1IBBc.19G>A (Hb C)SNV/indel2IBBc.20A>T (Hb E)SNV/indel5IBBc.20A>T (Hb S)SNV/indel3IBBc.126_129delSNV/indel3IBBc.126_129delSNV/indel2IBBc.137C>GSNV/indel1IBBc.316-197C>TSNV/indel1IBBc.316-197C>TSNV/indel1IBBc.82G>T (Hb Knossos)SNV/indel1	HBA	4.2del	HBA2 del	2	100%
HBA1 (Hb G-Philadelphia) 1 HBB (Hb Lepore-Baltimore) Exon 1, 2 del 2 100% HBB (Sicilian [deltabeta]0-Thal) Exon 1-3 del 1 100% HBB c.19G>A (Hb C) SNV/indel 2 100% HBB c.79G>A (Hb E) SNV/indel 5 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 2 100% HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c137C>G SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBA2		SNV/indel	1	
HBB (Sicilian [deltabeta]0-Thal) Exon 1-3 del 1 100% HBB c.19G>A (Hb C) SNV/indel 2 100% HBB c.79G>A (Hb E) SNV/indel 5 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 2 100% HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c.93-21G>A SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1	HBA1		SNV/indel	1	
HBB c.19G>A (Hb C) SNV/indel 2 100% HBB c.79G>A (Hb E) SNV/indel 5 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.20A>T (Hb S) SNV/indel 2 100% HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c137C>G SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	(Hb Lepore-Baltimore)	Exon 1, 2 del	2	100%
HBB c.79G>A (Hb E) SNV/indel 5 100% HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c.93-21G>A SNV/indel 1 100%	HBB	(Sicilian [deltabeta]0-Thal)	Exon 1-3 del	1	100%
HBB c.20A>T (Hb S) SNV/indel 3 100% HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c137C>G SNV/indel 1 100% HBB c138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	c.19G>A (Hb C)	SNV/indel	2	100%
HBB c.126_129del SNV/indel 2 100% HBB c.93-21G>A SNV/indel 2 100% HBB c137C>G SNV/indel 1 100% HBB c138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	c.79G>A (Hb E)	SNV/indel	5	100%
HBB c.93-21G>A SNV/indel 2 100% HBB c137C>G SNV/indel 1 100% HBB c138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	c.20A>T (Hb S)	SNV/indel	3	100%
HBB c137C>G SNV/indel 1 100% HBB c138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.32G>T (Hb Knossos) SNV/indel 1 100%	HBB	c.126_129del	SNV/indel	2	100%
HBB c138C>T SNV/indel 1 100% HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	c.93-21G>A	SNV/indel	2	100%
HBB c.316-197C>T SNV/indel 1 100% HBB c.82G>T (Hb Knossos) SNV/indel 1 100%	HBB	c137C>G	SNV/indel	1	100%
HBBc.82G>T (Hb Knossos)SNV/indel1100%	HBB	c138C>T	SNV/indel	1	100%
	HBB	c.316-197C>T	SNV/indel	1	100%
	HBB	c.82G>T (Hb Knossos)	SNV/indel	1	100%
	HBB	c.208G>A	SNV/indel	1	100%





ormal Intermediate Premutation




Figure 5. Predicted Risk of FMR1 Expansion Based on AGG Interruption Status. A cohort of 26 intermediate and premutation alleles were assessed using Asuragen Xpansion Interpreter[®] (XI) and PCR/nanopore using a custom algorithm. Genotypes were in 100% agreement with XI for the absolute number of AGG interruptions and within ±1 for the absolute position of each AGG interruption within the CGG repeat. Asterisks denote samples where AGG status modified the risk for a full mutation in the next generation compared to CGG repeat information alone.

 Table 4. Fourteen Carriers (5.7%) were Identified in a Presumed Normal Cohort
 (n=247) Using the PCR/nanopore Panel. CFTR T5011 is pending confirmation. FMR1 intermediate expansions were identified in seven samples. No HBB carriers were identified.

Sample ID	Carrier	FMR1	CFTR	SMN1, SMN2	HBA1/2
SID489	CFTR	30	R117H, WT	2,1	aa/aa
SID502	CFTR	23,30	G622D, WT	2,1	aa/aa
SID526	CFTR	20,29	T501I, WT	3,2	aa/aa
SID481	CFTR	31	F508del, WT	2,1	aa/aa
SID530	CFTR	30	F508del, WT	2,1	aa/aa
SID562	CFTR	30	F508del, WT	2,2	aa/aa
SID569	CFTR	30	F508del, WT	2,2	aa/aa
SID637	CFTR	30	F508del, WT	2,3	aa/aa
SID461	SMN1	20,30	WT, WT	1,1	aa/aa
SID546	SMN1	29	WT, WT	1,2	aa/aa
SID403	HBA2	30	WT, WT	2,1	3.7del/aa
SID415	HBA2	20,29	WT, WT	2,1	3.7del/aa
SID531	HBA2	29	WT, WT	2,2	3.7del/aa
SID557	HBA2	30	WT, WT	2,2	3.7del/aa

Genes for alpha-thalassemia, SMA, and FXS (HBA1/2, SMN1, and FMR1, respectively) are technically challenging due to homologous or GC-rich sequences. Repetitive polymorphisms, exon del/dups and other CF variants within CFTR can also be problematic to resolve. Here, we combine scalable PCR enrichment with nanopore sequencing in a prototype assay that accurately identifies multiple classes of pathogenic mutations for each of the seven carrier screening genes. The panel interrogates single nucleotide variants (SNVs), insertion/deletions (INDELs), exon del/dups, structural variants (SVs) and enumerates copy number variation (CNVs) and short tandem repeats (STRs) in a single workflow.

METHODS

Cell-line (CL) genomic DNA (gDNA) samples (n=168) were obtained from Coriell Cell Repository. Genomic DNA was also isolated from whole blood (WB) donors (n=249). Samples 3A were PCR amplified across three reactions, barcoded per sample, pooled across samples, and prepared using a ligation sequencing kit (LSK110 & LSK114; Oxford Nanopore Technologies, ONT). Sequencing was performed using MinION flow cells (R9.4.1, R10.4.1) on a Mk1B (ONT). Cellline samples representing all major classes of variants were used to develop custom data analysis pipelines and analyses were performed using custom software. Performance was demonstrated across cell-line and whole blood samples. Orthogonal methods or reporting (e.g. Coriell, 1000 Genomes, curve PCR analysis, custom PCR/capillary melt electrophoresis (CE), AmplideX[®] PCR/CE *CFTR* Kit[†], PCR/CE SMN1/2 Plus Kit⁺, and PCR/CE FMR1 Kit⁺, Xpansion Interpreter^{®‡} (XI), Sanger sequencing, and qPCR) were utilized to determine performance.

Figure 2. Sequence Data Differentiates SMN1 and SMN2 and Informs Silent Carrier Status and Disease Severity. A) Differentiation and assignment of reads to SMN1 and SMN2 by three different paralog-specific variants align reads to each gene and inform CNV prediction model. B-D) Silent carrier (SC1 or SC2) or disease modifier (DM) variant alignments.

SC2; c.*211_212del

g a c t c t a t t t t g a a

· · · · · · · ·

SMN1

Exon8




Figure 4. HBA Deletion Identification In Samples Based on Amplicon Coverage Using a Machine Learning Algorithm. Normalized fold change in coverage compared to control sample reference amplicons is shown on the y-axis; assay amplicons are arranged in sequential order on the x-axis. A) Typical WB sample with coverage pattern for two wild-type alleles (aa/aa). B) Typical WB silent carrier with one wild-type allele (aa) and one 3.7del allele leading to an α^+ -thal silent carrier where one copy of *HBA2* is effectively deleted (hybrid *HBA1/2* gene). **C)** CL sample with only one copy of *HBA1* and *HBA2* leading to a α^0 -thal trait. **D)** CL sample with an SEA deletion on one allele and an anti-3.7 triplication on the other. HBA2 and Rgn12 amplicons receive "wild-type" copy number calls because the combination of variants results in the same coverage levels as a wild-type sample (1 copy of HBA1, 2 copies of *HBA2*.


Table 2. CFTR Sample Level Agreement with Orthogonal Data for 103 Cell-Line and 234 Whole Blood Samples. The assay used Clair3² (SNV/indel) and read depth heuristics (del/dup) to detect 57 unique variants among 97 total variants, including two del/dup (CFTRdele2,3, CFTRdele19-21), which represent 88.9% prevalence of variants in the ethnically diverse US population.¹ The assay was performed with both R10.4.1 (n=290), and R9.4.1, (n=47). [¥]Sanger sequencing verification pending for two variants (CL, c.3368-2A>T; WB, T501I), with corresponding nanopore read pile-ups shown on the right.

Cell Line		Orthogo	nal Data		Overall	TCTTTGGGGATCTAT HET WT	10 bp 117,814,828 bp 1
Sample Level Accuracy	Allele1/Allele2	WT/WT	MUT/WT	MUT/MUT	Genotype Agreement		:
pore	WT/WT	48	0	0	50/50 (100%)) (c.3368	_
<i>CFTR</i> /nano	MUT/WT	0	27	0	27/27 (100%)	A12960	· · · ·

CONCLUSIONS

- The prototype PCR/nanopore assay accurately resolves multiple challenging variants across several classes for seven of the most common carrier screening genes.
- At least 96 barcodes can be combined in a single run using three PCR enrichment reactions and a unified data processing pipeline.
- The assay utilizes a single-platform, streamlined workflow and has potential to greatly reduce carrier screening complexity and turn around times, especially when paired with other similar assays (see Poster #P455, P546).
- In over 400 samples tested across each gene, the PCR/Nanopore assay agreed with the orthogonal method for SNVs/INDELs in SMN1, CFTR, HBA1, HBA2, and HBB (>99% of samples), SMN1 CN (96.6%), SMN2 CN (98.5%), FMR1 repeat categories (100%), FMR1 AGG interruptions (100%), and HBA1/2 SVs (100%).

Figure 3. SMN1/2 PCR/Nanopore Assay Classifies Carrier with 100% accuracy. Calling accuracy for SMN1 and SMN2 copy numbers in A) 96 CL and B) 227 WB samples. Hyperparameters for the decision tree model were selected using an 80:20 train:test split in a stratified randomly selected five-fold cross validation scheme on an independent set of 349 samples (62 CL and 287 WB). All carriers were identified with 100% accuracy (green dashed outline).



REFERENCES

Beauchamp, K. A. et al. (2019) Genet Med 37, 773-8.

2. Zheng, Z. et al. (2022). Nat Comp Sci 2, 797–803.



This product is under development. Future availability and performance to be determined.

[†]For Research Use Only. Not for use in diagnostic procedures. [‡]Xpansion Interpreter[®] is a laboratory-developed test. Analytical and clinical performance have not been reviewed by the FDA. All authors have the financial relationship to disclose: Employment by Asuragen. Presented at ACMG 2023

2B

SC1; c.*3+80T>G

T T T T G A A C A T T T A A A A A A

* * * * *

SMN1

Intron7