Accurate Genotyping of Complex, Clinically-Relevant Variants in 11 Hard-to-Decipher Genes by **Combining Novel PCR Chemistries with Any-length Nanopore Sequencing for Carrier Screening**

Bradley Hall, Bryan Killinger, Christopher J. Fraher, Cody Edwards, Ryan Routsong, Jonathan Turner, Monica Roberts, Jon Kemppainen, Brennan Greenlee, and Gary J. Latham Asuragen, a Bio-Techne Brand, Austin TX

Summary

- Cystic Fibrosis (CF), Spinal Muscular Atrophy (SMA), Hemoglobinopathies, and Fragile X Syndrome (FXS) are among the most inherited genetic disorders.
- High-prevalence carrier genes, associated with disorders such as Gaucher Disease (GD), congenital adrenal hyperplasia (21-OHD CAH), and Hemophilia A (HA), include complex structural variants and pseudogenes that confound conventional sequencing methods.
- We explored whether novel PCR enrichment, nanopore sequencing, and machine learning models could detect multiple classes of variants including SNVs, INDELs, Exon del/dups, SVs, CNVs and STRs in a single workflow.
- The assay was optimized with 233 cell-line samples and performance evaluated with 347 whole blood, 1 saliva, and 5 buccal samples across the eleven genes with a subset set of 233 whole blood to identify potential carriers from presumed normal donors.

Introduction

Nearly everyone is a genetic carrier of a disease or condition¹. Carrier screening (CS) identifies couples at risk for having a child with a severe genetic disorder. Although Next-Generation Sequencing (NGS) is a widely used method, it fails to resolve many problematic genes recommended in professional practice guidelines due to GC-rich tandem repeats, copy number variation, pseudogenes, and structural variation. More broadly, 20.4% of pathogenic/likely pathogenic variants in ClinVar have been reported in "dark" or "camouflaged" regions of the genome that are "technically challenging" to resolve². Many of these genes require specialized techniques and only cover a fraction of carrier risk.

To address these shortcomings, we combined three innovations : 1) short- & long-range PCR, 2) any-length nanopore sequencing, and 3) customized software analysis pipelines. Using a single workstream, we developed a modular panel of 11 genes critical for CS, including nine "hard-todecipher" genes. We estimate the panel represents ~75% of at-risk couple detection compared to panels that are 15 times larger³. Here we describe results utilizing this prototype assay to genotype CFTR, SMN1, SMN2, FMR1, HBA1, HBA2, HBB, F8 intron inversions, GBA, CYP21A2, and TNXB across more than 400 samples.

Methods

Cell-line (CL) genomic DNA (gDNA) samples (n=233) were obtained from Coriell Cell Repository. Genomic DNA was also isolated from presumed normal or clinical donors (n=353). One contrived and two commercial control samples were utilized for F8 intron inversion. Target regions were enriched in four PCR reactions, barcoded per sample, pooled across samples, and prepared using a ligation sequencing kit (LSK110 & LSK114; Oxford Nanopore Technologies, ONT). Sequencing was performed using MinION flow cells (R9.4.1, R10.4.1) on a Mk1B (ONT). A cohort of 12 cell lines and 232 presumed normal whole blood samples were analyzed by the entire panel. Remaining samples were tested by a subset of primer mixes or genes analyzed for variants. Cell-line samples representing all major classes of variants were used to develop custom data analysis pipelines and software. Clair3 was utilized for SNV/INDEL identification⁴. Performance was demonstrated across cell-line and whole blood samples. Orthogonal methods or reporting (e.g. Coriell, 1000 Genomes, melt curve PCR analysis, MLPA, custom PCR/capillary electrophoresis (CE), AmplideX[®] PCR/CE CFTR Kit[†], PCR/CE SMN1/2 Plus Kit[†], and PCR/CE FMR1 Kit[†], Xpansion Interpreter^{®‡} (XI), Sanger sequencing, and gPCR) were utilized to determine comparator results.



Figure 1. PCR/nanopore Carrier Screening Panel Design and Workflow Identifies Pathogenic Variants for 11 Genes Responsible for Common Inherited Genetic Disorders and/or Challenging to Assay by Convention Methods. The combination of AmplideX[®] PCR technology across **A**) 4 primer mixes, and **B**) nanopore sequencing enables detection of **C**) multiple variant classes for each of the 11 genes within the panel in a single workflow.

Results



Differentiation and assignment of reads to SMN1 and SMN2 by three different paralog-specific variants align reads to each gene and inform CNV prediction model. B) Silent carrier (SC1; c.*3+807>G) variant alignments. SC2 (c.*211 212del) and disease modifier (DM; c.859G>C) not shown.



Figure 3. SMN1/2 PCR/Nanopore Assay Accurately Classifies Carrier Status. Calling accuracy for SMN1 and SMN2 copy numbers in A) 97 CL and B) 232 presumed normal and 70 clinical WB samples. Hyperparameters for the decision tree model were selected using an 80:20 train:test split in a stratified randomly selected five-fold cross validation scheme on an independent set of 349 samples (62 CL and 287 WB). C) Carriers were identified with 100% accuracy (green dashed outline). There were 3 clinical samples removed from analysis for *SMN2* due to lack of orthogonal data.



Figure 4. CFTR Sample Level Agreement with Orthogonal Data for 104 Cell-Line and 234 Whole Blood Samples. The assay used Clair3 (SNV/indel) and read depth heuristics (del/dup) to detect 65 unique variants, including three del/dup (CFTRdele2,3, CFTRdele19-21, CFTR dele4-11). A-B) The assay was performed with both R10.4.1 (n=291), and R9.4.1, (n=47) with 100% agreement or orthogonal methods. **C-D)** Nanopore read pile-ups of two variants previously undetected in whole blood (G622D, T501l; [¥]Sanger sequencing verification pending).



Figure 5. Large Exon Deletions Detected at 100% PPV and Sensitivity. Normalized amplicon coverage distinguishes exon deletion (orange) from wildtype genotypes (teal) for A) CFTR2,3dele and B) CFTRdele20. Samples with normalized coverage below the threshold (vertical dashed line) were classified as heterozygous for the large exon deletions.

Table 1. Pathogenic and Likely Pathogenic *HBA1/2* and *HBB* Variants Detected in 147 Cell-Line and 283 Whole Blood or Clinical Samples. The sample set contained 247 presumed normal WB, 30 clinical WB, 5 buccal, 1 saliva tested with the combined HBA1/2 and HBB assay on R10.4.1 flow cells. Of these, 94 CL and 238 WB samples were genotyped as wild-type (aa/aa). All but 3 samples were concordant with orthogonal data. The three discordant samples were miscalled for anti-3.7 (2 FN, 1FP) and do not affect carrier status.

Gene	Variant Type	Unique Variants Detected	Sample Agreement
HBA	Structural	SEA, FIL, THAI, 3.7del, anti-3.7, and 4.2del	72/75 (96%)
	SNV/Indel	c.427T>C (Hb Constant Spring), c.207C>G (Hb G-Philadelphia), c.237del, c.142G>C, c.95+2_95+6del	8/8 (100%)
HBB	Structural	(Hb Lepore-Baltimore), (Sicilian [$\Delta\beta$ 0-Thal), HPFH-1, HPFH-2	8/8 (100%)
	SNV/Indel	c.19G>A (Hb C), c.79G>A (Hb E), c.20A>T (Hb S), c.126_129del, c.93-21G>A, c137C>G, c138C>T, c138C>T, c.118C>T, c.315+1G>A, c.316-197C>T, c.82G>T (Hb Knossos), c.316-2A>G c.17_18del, c79A>G, c.404T>C, c.27dup, c.364G>C, c.92+5G>C, c.92+6T>C, c.92G>C, and c.208G>A	44/44 (100%)



Predicted

Table 2. FMR1 Categorical Agreement with Orthogonal Genotypes for 99 Cell-Line and 232 Whole Blood Samples. Using R9.4.1, ACMG categorical genotype boundaries are included for reference. All samples fell within expected categories based AmplideX PCR/CE FMR1 precision metrics (±1:0-70 repeats, ±3:71-119). All expanded samples, including full mutations up to 940 CGG repeats, were flagged correctly. Additionally, CGG sizing was accurate within precision for 321/331 (97.0%) samples and 443/454 (97.5%) alleles. In 7/11 samples, the algorithm accurately called one of two alleles when two similar sized alleles (1-3 repeats difference) were present. The algorithm identified a previously unidentified minor mosaic allele in the remaining three samples that did not change the categorical call.

Sample level categorical accuracy	Normal <45 CGG	Intermediate 45-54 CGG	Premutation 55-200 CGG	Full Mutation >200 CGG	Sensitivit
Training	61	18	55	19	100%
Cell line	76	5	14	4	100%
Whole Blood	225	7	0	0	100%



Figure 6. Predicted Risk of FMR1 Expansion Based on AGG Interruption Status. A cohort of 26 intermediate and premutation alleles were assessed using Asuragen Xpansion Interpreter[®] (XI) and PCR/nanopore using a custom algorithm. Genotypes were in 100% agreement with XI for the absolute number of AGG interruptions and within ±1 for the absolute position of each AGG interruption within the CGG repeat. The cyan dots denote samples where AGG status modified the risk for a full mutation in the next generation compared to CGG repeat information alone.



Figure 7. Accurate Resolution of Copy Number and Pseudogene Fusions in the CYP21A2 Gene Cluster Utilizing Sequencing **Deconvolution on R10.4.1**. 28 pathogenic variants were identified across cell-line and whole blood samples and mapped to their location in *CYP21A2*. We did not find evidence of *TNXB* in samples tested.

Figure 8. Prototype Assay Amplifies and Identifies >100 kb F8 Inversions. F8 inversions were observed in three samples using R10.4.1. These include a homozygous intron 22 inversion "affected", a heterozygous intron 22 inversion "carrier", and a contrived heterozygous intron 1 inversion "carrier" created by mixing an inversion g-Block and a wild-type sample. The results show the proportion (y-axis) of *F8* reads that aligned to either the reference genome (Wild Type) or an F8 inversion. Inversion reads were not identified in the wild type sample. Both intron 1 and intron 22 inversions were verified by orthogonal methods (data not shown).

Figure 9. Pathogenic and Likely Pathogenic HBA1/2 and HBB Variants Detected Using the PCR/nanopore Assay. HBA deletion identification was based on normalized fold change in read depth compared to control sample reference amplicons on the y-axis; assay amplicons are arranged in sequential order on the x-axis. WB sample with a FIL deletion on one allele

10A. GBA_Ex01-11_T21408_T21409_0 BA_Ex01-11_T21408_T21409_1



Figure 10. Sequence Data Reveals Diverse Pathogenic Variants Across the GBA gene on R10.4.1. Clinically affected cell-line NA20270 is a compound heterozygote and contains two SNVs within the GBA gene; A) one allele has a T>C transition at nucleotide 1448 in exon 10 (L444P. c.1448T>C), and B) the other allele has a splice site mutation in intron 2 (IVS2+1G>A). The variants in the red box are pseudogene microconversions consistent with GBAP1 paralog-specific variants. All samples were verified with orthogonal methods (data not shown).









10B.









Figure 11. PCR/nanopore Carrier Screening Panel Design and Workflow Identifies Pathogenic Variants for 11 Genes Responsible for Common Inherited Genetic Disorders and/or Challenging to Assay by Convention Methods. The combination of AmplideX® PCR technology across 4 primer mixes, and nanopore sequencing enables detection of multiple variant classes for each of the 11 genes within the panel. A total of 589 samples were utilized for training and testing. For a subset of samples, only specific gene data was analyzed and compared, especially if a variant was known in the gene of interest.

Table 3. Detected 37 Carriers (15.9%) in a Presumed Normal Cohort of Whole Blood Samples (n=232) using the Prototype **Assay.** One donor sample was identified as a carrier for both *SMN*1 and *CYP21A2*. *CFTR* T5011 is pending confirmation. *FMR*1 intermediate expansions were identified in seven samples but not shown. Hemophilia A (F8 intron inversion) carriers were not identified (expected by carrier rate). All CYP21A2 and GBA variants are pathogenic or likely pathogenic, however the lighter grey CYP21A2 variants are non-classical and often go under-diagnosed. All variants were confirmed by orthogonal methods. *SC variants suggest an increased carrier risk.

Gene	# of Carriers	Variants Identified
CFTR	8	5x F508del, wt; R117H, wt; G622D, wt; T501l, wt
SMN1	2	1, 1 and 1, 2 (<i>SMN</i> 1, <i>SMN</i> 2)
SMN1, SC	3	2, 2 + SC1, SC2 [†]
HBA1/2	4	3.7del/aa
HBB	3	c.404T>C; c.20A>T (Hb S); c79A>G
GBA	3	2x N370S (c.1226A>G); R496H (c.1604G>A)
CYP21A2	4 8	30-kb del; 3x Q318X (c.955C>T); 7x V281L (c.844G>T); P453S (c.1360C>T)
SMN1 & CYP21A2	1	1,3 & Q318X (c.955C>T)

Conclusions

- The prototype PCR/nanopore assay accurately resolves genotypes across multiple classes of variation for 11 of the most common and challenging genes associated with heritable disease.
- The assay utilizes a single-platform, streamlined workflow, and has potential to greatly reduce carrier screening complexity and turn around times.
- Detection of a dual carrier (CYP21A2 and SMN1) demonstrates the potential of a unified workflow to address targets that are conventionally analyzed through disparate assays.
- In 589 samples tested across the panel, the PCR/nanopore assay agreed with the orthogonal methods for SNVs/INDELs in SMN1, CFTR, GBA, CYP21A2, HBA1, HBA2, and HBB (>99% of samples), SMN1 CN (97.2%), SMN2 CN (98.2%), FMR1 repeat categories (100%), FMR1 AGG interruptions (100%), and *HBA1/2* deletions (100%).

We would like to thank Jaime L. Lopes, Marcus Cannon, Charles T. Quinn, and D. Brian Dawson from Cincinnati Children's for providing residual clinical samples and orthogonal method results (HBA/HBB). We would like to thank Ahmad Tayoun and Sawsan Yaslam from Al Jalila Children's Specialty Hospital, Dubai, UAE for testing residual clinical samples and orthogonal method results (SMN1/2).

[†]This product is under development. Future availability and performance to be determined. [‡]Xpansion Interpreter is a laboratory-developed test. All authors have the financial relationship to disclose: Employment by Asuragen.

Presented at AMP 2023

References:

- 1. Xiao, Q. et al. Npj Genom Medicine 6, 41 (2021).
- 2. Lincoln SE., *et al*. Genet Med 23 (2021).
- 3. Beauchamp, KA., *et al*. Genet Med 21 (2019).
- 4. Zheng, Z., *et al*. Nat Comp Sci 2, 797–803 (2022). 5. Taber, K. J. *et al*. Genet Med (2021).



G030